



Ompi
Stevanato Group

**How to assure quality
in glass vials:
controlled
manufacturing
processes focused on
particle reduction**

Martina Largoni
Product Manager Vial Platform



How to assure quality in glass vials

- **Stevanato Group at a glance**
- Bulk manufacturing process optimization
- State-of-the-art solutions for EZ-fill products
- Visible Particle Reduction Program
- Conclusions

Stevanato Group Brand Structure



PHARMACEUTICAL SYSTEMS

ENGINEERING SYSTEMS


SERVICES

Stevanato Group Today – Global Footprint



Ompi Range of Products

Bulk Containers




.... A Complete Range of Containers for Injectables

A blue decorative graphic in the top right corner, consisting of a horizontal bar with a diagonal cutout and a vertical extension on the right side.


How to assure quality in glass vials

- Stevanato Group at a glance
- **Bulk manufacturing process optimization**
- State-of-the-art solutions for EZ-fill products
- Visible Particle Reduction Program
- Conclusions

Glass Forming Process: from Bulk to EZ-fill®



Bulk Products | Technology Steps and Improved Solutions



Reduction of Variability in Dimensions

CONTINUOUS IMPROVEMENT OF THE TECHNOLOGY



Forming tools are designed to reduce glass container tolerances and to maintain their precision for long forming runs

INCREASED NUMBER OF FORMING STEPS



Forming steps are designed to guarantee high precision forming and high repeatability of the process

NEW GENERATION INSPECTION SYSTEM



100% in-line camera inspection gives a real-time feedback on the quality of the batch



Example of Improved Manufacturing Process in Vial Bottom Forming Step

- Increased number of forming steps
- Introduction of a mold for final manufacturing
- Pyrometer Technology for 100% temperature control


Higher stability of vials

Homogeneous glass distribution
and small concavity

More resistant glass in
lyophilization cycles



Accurate Handling to Preserve Cosmetic and Mechanical Properties of Glass Container



- No glass-to glass contact
- No buffer stations
- New contact materials to reduce the risk of thermal shock and avoid metal to glass
- Soft handling of the glass container to limit vibrations

100% in-line Inspection Controls are in Place to Assure Product Conformity

- 100% inspection of all dimensions
- Automatic rejection of defective pcs
- Automatic calibration system
- Performance per chuck
- Measurements and statistics in real time




Specific Line Settings Contribute to Achieve Different Quality Levels

	SG Ompi Fina	SG Ompi Nexa
Critical		
Chip (sealing) Crack Glass particles	0,1 – 0,04	0,1 – 0,025 & ppm
Major		
Contamination Bull eye Groves and notches	0,65 – 0,1	0,15 – 0,04
Minor		
Notch, external Scratch Folds	2,5 – 0,4	1,0 – 0,25

Glass Chemical Properties

Glass-Liquid Interactions Can Lead to Creation of Altered Layer



Alkaline solutions strongly affect the dissolution of the silica layer.
 SiO_2 concentration in the extraction liquid increases steeply

Flakes appears

Several Factors Affect Delamination Propensity of Pharmaceutical Glass



Morphological and Physicochemical properties of Glass Vials can affect the interaction with the drug

Container manufacturer direct influence

Glass Vial

Conversion manufacturing process
(e.g.: speed, temperature, type of glass)

Surface treatments
(e.g.: sulfur)

Pharma company direct influence

Drug

Formulation
(e.g.: chemical composition, pH)

Process
(e.g.: depyro, sterilization)

Storage conditions
(e.g.: temperature, time)

Source: USP 1660

Chemical Performances Can Be Guaranteed with Optimized Bulk Processes

LDP (Low Delamination Propensity) Vials



Forming process optimization with low heat/energy thermal cycle and reduced surface inhomogeneities formation



Quantitative and qualitative **tests** to guarantee the quality and the stability of vials production




No coatings
No glass formulation changes
No need to re-file

Responsibilities

- Primary Packaging Supplier optimizes glass Converting Process
- Pharma Companies in charge of verifying impact of washing/depyrogenation

LDP Vials Show Chemical Performances Close to Non-Converted Raw Material



Example of extracted elements from converted glass in comparison to raw materials


Vials from Optimized Process shows comparable performances of Raw Material
Converting process can be controlled to reduce effect on chemical properties of the glass




How to assure quality in glass vials

- Stevanato Group at a glance
- Bulk manufacturing process optimization
- **State-of-the-art solutions for EZ-fill products**
- Visible Particle Reduction Program
- Conclusions

Ompi EZ-fill® Syringes, Vials and Cartridges | Concept Introduction



OMPI Designed Different Secondary Packaging Configurations to Support the Specific Needs of Pharma



↓

Washing
(Siliconization)
Depyrogenation
Nesting (no G2G)
Final Sterilization



OMPI EZ-fill in **Tray** is a suitable solution for EZ-fill vials and cartridges for:

- High speed filling line
- Flexible lines processing different dimensions of the same container
- RTU lines working according to traditional filling process (i.e. for cartridges: plunger insertion, filling from the neck, capping)



OMPI EZ-fill in **N&T** is the best configuration for *combi lines* to fill vials, cartridges, syringes with a unique filling machine

Ompi EZ-fill® | Pharma Production Area

Production Area

Designed and developed on GMP standard.



Different Clean Rooms

(from ISO8 to ISO5)




Internal Laboratories

(chemical, environmental and functional testing)



EZ-fill® Products | Technology Steps and Microbiological Focus



Product Microbiological Level

- Endotoxin Level < 0,25 EU/mL
- Bioburden (before sterilization)
- Final Sterility < 0 CFU

Utilities Microbiological Periodical Analysis

- Water for Injection
- Compressed Air
- Environmental monitoring (e.g. particle counting)



How to assure quality in glass vials

- Stevanato Group at a glance
- Bulk manufacturing process optimization
- State-of-the-art solutions for EZ-fill® products
- **Visible Particle Reduction Program**
- Conclusions

Visible Particle Characterization on 10R EZ-fill® Vials in Nest&Tub

Specification on primary packaging prescribes a requirement for visible particles bigger than 0.3 mm^2

Focus on particles/fibers **bigger than $100 \mu\text{m}$**

Characterization results on 2 tubs of 10R vials manufactured in 2016:

- Identified fibers (L>3W): Cellulose, Polyester.
- Identified particles: Proteins, Cellulose, Teflon, Poly-Acetal.
- No secondary-packaging-related materials were identified (e.g. Polypropylene, Tyvek, Polyethylene).



Samples were characterized through automatic scanning on filter (USP <788>), using a rinsing solution 0.2 μm -filtered water (considered particle-free), added with 0.001% Tween 80.

Critical-to-Quality Attribute (CTQ's)

1. CTQ (OPERATIVE): Defect per Unit (DPU)

$$CTQ_{ope} := DPU = \frac{N_{VP,tot}}{N_{vials,inspected}}$$

where:

- DPU is the mean number of visible particles per vial;
- $N_{VP,tot}$ is the total number of «generalized» particles (as sum of particles and fibers);
- $N_{vials,inspected}$ is the total number of inspected vials.

2. CTQ (PROJECT): Process Yield (PY)

$$CTQ_{prj} := PY = e^{-DPU}$$


where e is the Euler's number.

ASSESSMENT BEFORE IMPROVEMENT			
Particles/ Fibers size > 100 µm (Study threshold)		Particles/ Fibers size > 0.3 mm ² (Specification threshold)	
DPU	PY [%]	DPU	PY [%]
>1	<50%	0	100.00%


Test plan: two (2) random tubs from 2016 batch (48 vials/tub) – tot. 96 vials

Manufacturing Flow Optimization

BEFORE IMPROVEMENT



AFTER IMPROVEMENT



- Dedicated semifinished pallet loading area
- Improved material movimentation
- Additional changing rooms

Performance Improvements Based on Manufacturing Flow Optimization: Effectiveness Verification

ASSESSMENT BEFORE IMPROVEMENT			
Particles/ Fibers size > 100 μm (Study threshold)		Particles/ Fibers size > 0.3 mm^2 (Specification threshold)	
DPU	PY [%]	DPU	PY [%]
>1	<50%	0	100.00%

Test plan: two (2) random tubs from 2016 batch (48 vials/tub) – tot. 96 vials



SG OMPI PROCESS MAPPING			
Particles/ Fibers size > 100 μm (Study threshold)		Particles/ Fibers size > 0.3 mm^2 (Specification threshold)	
DPU	PY [%]	DPU	PY [%]
<0.03	>97%	0	100.00%

Test plan: twelve (12) tubs from 2017 production process mapping (48 vials/tub) – tot. 12 tubs (576 vials)



How to assure quality in glass vials

- Stevanato Group at a glance
- Bulk manufacturing process optimization
- State-of-the-art solutions for EZ-fill® products
- Visible Particle Reduction Program
- **Conclusions**



Conclusions

- High quality level for glass primary packaging can be obtained starting from the optimization of manufacturing process
- The characteristics obtained during bulk forming process can be maintained through the EZ-fill step thanks to a state-of-the-art design of the lines
- Particle Reduction is possible acting on flow optimization in core area
 - Further improvements are possible through additional assessment on raw material Suppliers' processes
- The Ready-to-Use concept is today well established in the market since it allows to:
 - Reduce the total cost of ownership
 - Increase flexibility
 - Increase quality



Thank You for Your Attention!

For further information visit
www.sg-ompi.com