PROCESSING & PACKAGING

RABS vs. Isolator – Status Quo and a look into the future

M. Angelmaier, Global Product Manager Isolator Technology

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3

Meaning

Green emphasizes the importance of sustainability and health.

Technology is our core competence.

The square symbolizes product protection through packaging technology.

SYNTEGON

Synergy

Today's challenges are too big to solve alone. We rely on working in partnership with our customers, suppliers and experts.

On

We look forward and act in a future-oriented way.

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Agenda

5

- 01 Introduction
- 02 Characteristics of RABS and Isolator systems
- 03 Current and future developments
- 04 Summary and Outlook

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
Introduction

6

Introduction – Aseptic Risk Reduction

Barrier systems like RABS or Isolators offer “tangible advantages over traditional aseptic processing” - FDA, *Aseptic Guidance, Appendix 1*

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
Attributes and Layouts – Open RABS


“Active” RABS

- Integrated UDAF with diffuser screen or plate
- HEPA filters for air filtration

“Passive” RABS

- Room UDAF
- Air filtration & homogenization through clean room ceiling

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Introduction


9

Attributes and Layouts – Closed RABS


“Active” RABS

- Integrated AHU with diffuser screen or plate
- HEPA filtration / Air Recirculation

“Passive” RABS

- Room AHU
- Air recirculation
- Filtration by clean room ceiling

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
SYNTEGON

Attributes and Layouts – Isolator

- **Isolator design features:**
 - (Hard)Wall between operator and aseptic area
 - Positive pressure “inside vs. outside”
 - Sterile transfer of product, container and components
 - HEPA14 filtered air
 - Uni-Directional Air Flow
 - Reproducible decontamination (system)
- Class C or D surrounding cleanroom

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Regulatory

RABS

In the past RABS or closed RABS are not clearly defined in any standard, guidance or regulation.

- However, a lot of requirements – particularly for maintaining ISO class 5 above open container – are also mandatory for RABS.
- In 2005 FDA asked ISPE to create a definition of 'RABS'. The paper was developed within a team of industry experts together with FDA (Rick Friedman).
- The PHSS has worked out a Technical Monograph on RABS which is accompanied by a RABS classification.

Isolator

Clearly defined regulations exist

- EU GMP Guideline for „Medicinal Products“, Annex 1
- ISO 13408 „Aseptic Processing of Health Care Products“, Part 6
- EN 13824 „Sterilization of Medical Devices“
- PDA Guideline „Design and Validation of Isolator Systems...“
- ISPE Baseline Guide™ „Sterile Manufacturing Facilities“
- PIC/S „Recommendations on the Inspection of Isolator Technology“
- FDA Guidance for Industry „Sterile Drug Products Produced by Aseptic Processing“
- DIN EN ISO 14644-7 „Cleanrooms and associated controlled environments“, Part 7: Separative Devices

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SYNTEGON

Regulatory

- Within current Annex1 draft (February 2020) both **RABS** as well as **Isolator** systems are mentioned
- Traditional cleanroom installations are seen critical for new installations „...Any alternative approaches to the use of RABS or isolators should be justified...“


181 4.3 Restricted Access Barrier Systems (RABS) and isolators are beneficial in assuring the required
182 conditions and minimizing the microbial contamination associated with direct human interventions
183 in the critical zone. Their use should be considered in the CCS. Any alternative approaches to the use
184 of RABS or isolators should be justified.

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Cleaning and Decontamination

RABS

- Manual cleaning and disinfection
- Room biodec with H_2O_2 , NO_2 or other

Isolator

- Isolator biodec with integrated H_2O_2 system
- Usage of evaporated or nebulized H_2O_2

SafeVAP biodecontamination system

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Material Transfer


RABS

- ▶ Transfer airlocks for component supply
- ▶ Manual procedures

Isolator

- ▶ Fast H_2O_2 Transfer chamber
- ▶ Integrated RTP systems
- ▶ Usage of robotic systems for automated RTP handling

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
SYNTEGON

Energy Consumption

- **Huge difference for media consumption** depending on system and room conditions
- Syntegon Article in **Pharmaceutical Engineering** "Comparing Energy Consumption of RABS and Isolator Configurations"

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Pros/Cons Isolator vs. RABS

RABS

- Much lower capital expenditure
- Every day multi product manuf.
- Easier to adapt coming from CR
- Less qualification effort necessary
- Retrofit possible
- Lower SAL
- Not flexible for highly potent
- Operator can intervene
- High running costs

Isolator

- Highest Sterility Assurance Level
- Long campaigns
- Highly potent products
- Operator convenience
- Saves clean room space
- Lower running costs
- Highest capital invest
- Biodecontamination time
- Manual interventions difficult
- More complex process troubleshooting

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SYNTEGON

Current and future Developments – Usage of gloves

- Glove testing procedures also for RABS
- Usage of EPDM and CSM gloves
 - EPDM with additional white layer
→ improved visual glove inspection by operators
- Implementation of glove management systems

5.16 The design of the RABS or isolator shall take into account all critical factors associated with these technologies, including the quality of the air inside and the surrounding area, the materials and component transfer, the decontamination, disinfection or sterilization processes and the risk factors associated with the manufacturing operations and materials, and the operations conducted within the critical zone.

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

SYNTEGON

Current and future Developments – Material Transfer

- **More „robust“ procedures** for getting stoppers, caps, etc. into **Grade A area**
- Example for Integration of ATEC System for Stopper and Cap Transfer for RABS Application

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
SYNTEGON

Current and future Developments – Cleaning/Decontamination

- Integration of reproducible cleaning and decontamination methods for **RABS**
- Example „Vetter V-CRT“
- Usage of optimized H_2O_2 biodec for RABS and surrounding room

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SYNTEGON

Current and future Developments – Cleaning/Decontamination

- Optimization of biodec cycle times for **Isolator** systems and air locks
 - Total Cycle Time **down to 0,1 ppm in approx. 80 minutes** (1h 20min)
 - Total Cycle Time **down to 0,05 ppm in approx. 110 minutes** (1h 50min)
- “From hours down to minutes”
- Challenge for low residual H_2O_2 concentrations

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

SYNTEGON


Current and future Developments – CFD Simulations

- Usage of **CFD** simulations for flow-optimized machine designs and critical areas
- Implementation in early project phases

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

Current and future Developments – High-potent products

- Industry trend:** Increasing amount of oncological products in the pipeline
- Often high-potent API's with highest toxicity levels up to **OEB 6 with $<0,1 \mu\text{g}/\text{m}^3$**

Elements of an Antibody-Drug Conjugate (ADC)

References: 1. Carter PJ et al. *Cancer J*. 2008;14(3):154-169. 2. Senter PD. *Curr Opin Chem Biol*. 2009;13(3):235-244. 3. Polson AG et al. *Cancer Res*. 2009;69(8):2358-2364.

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30

Summary and Outlook

- Higher degree of automation and robotic systems
- “Gloveless” Barrier Systems
- Integration of robust and reproducible procedures for RABS:
 - Material Transfer
 - Cleaning and Biodecontamination
 - ...
- Optimized changeover for Isolators with rapid biodecontamination cycle times

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© Syntegon Gruppe. Alle Rechte vorbehalten, auch bzgl. jeder Vervielfältigung, Verbreitung, Reproduktion, Bearbeitung, Weitergabe sowie für den Fall von Schutzrechtsanmeldungen.