

Advancements on implementation of single use technology in vaccine manufacturing

October 25th, 2013

**DCVMN Meeting
Rio de Janeiro, Brazil**

Overview of single use products

Applications of single use technology in vaccine manufacturing

Validation of single use systems

Overview of SU Products

- SU products
 - Single use containers or bags for storage and sampling
 - Filter assemblies with tubing (silicone and /or C-Flex)
 - Non-sterile and sterile connectors
 - Disposable mixing system for solution preparation and formulation
 - Final filtration assemblies for final filling
 - Single use systems




Film Technology

Typical 2D Bag Subassemblies



- 50mL to 50L

Typical 3D Bag Subassemblies

- 100L to 3500L for any bins
- Hand free filling available



Film:



Sterile connectors

Description & functionality



What it is:

**Single-use sterile connector made with
1 female coupling + 1 male coupling part**

What it does:

An operator independent, sterile connection between

γ -sterilized assembly + γ -sterilized assembly

γ -sterilized assembly + Autoclaved assembly

Autoclaved assembly + Autoclaved assembly

... in any environment !



SU mixing systems and mix bags

Wide range of sizes (10/50/100/200/500/1000 L)

- Levitating, magnetically driven impeller
- Configurable bag assemblies

Electronic drive unit and motor

- Portable & removable
- Interchangeable
- Multiple carrier options
- PE, SS, Jacketed SS, Load Cells
- Stable & Mobile



Sampling system

Increased sampling productivity, while reducing set-up, cleaning and flushing time


Traditional Sampling Method

- Lengthy CIP/SIP with each sample
- Additional steps over and over again




Sampling with Millipore

- Only one SIP
- One to 45 samples at once
- Push to sample
- Unlimited possibilities
- Separate & crimp



Increasing Use of Single-use systems

Reasons for increasing Disposables in 2011
 (% Indicating Attribute is "Very Important" or "Important")




Overview of single use products

Applications of single use technology in vaccine manufacturing


Validation of single use systems

Single use technology in Vaccine Manufacturing



Media/Buffer Preparation


Use of disposable mixing systems with liners or bags.
Bag with filter assemblies
Media addition assemblies



Transfer lines and manifolds



Sterile sampling from SST bioreactors



REVISIONS			
REV.	DESCRIPTION	DATE	REVISED BY
1	NEW QUOTE	12-MAY-11	C. RAYE

ITEM NO.	DESCRIPTION	LG.	QTY.
1	CONN-LYNK ST 1/4"HB X 3/4"TC		1
2	BAG-NOVASEPTUM 250ML (5.12"W X 7.3"HT) END-PORTED (2 X 1/8")		6
3	BAG-NOVASEPTUM 1000ML (6"W X 12.3"HT) END-PORTED (2 X 1/8")		3
4	TUBING-PHARMA 65 (3mm ID X 6mm OD)		4
5	TUBING-PHARMA 65 (3mm ID X 6mm OD)		6
6	TUBING-CFLEX 374 (1/8"ID X 1/4"OD)		6
7	TUBING-CFLEX 374 (1/8"ID X 1/4"OD)		12
8	TUBING-CFLEX 374 (1/4"ID X 1/2"OD)		12
9	NOVASEAL FOR 3MMID X 6MMOD TUBING		8
10	COUPLER-CROSS 1/8"		4
11	COUPLER-STRAIGHT REDUCER 1/4" X 1/8"		1
12	CONN-LUER FEMALE 1/8"HB		9
13	CAP-LUER MALE KYNAR		9
14	CLAMP-PINCH 1/2" (6MM-1/8"OD TO 3/8"OD TUBE)		9
15	DUST COVER 6" X 6"		1
16	CERTIFICATE OF QUALITY - SILVER (NOT SHOWN)		3
17	LABEL - GAMMA STERILIZED ASSEMBLIES (NOT SHOWN)		3

NOTES:
1. ALL PARTS TO BE SUPPLIED AS A NON-STERILE UNIT BEFORE ASSEMBLY.
2. PLACE ASSEMBLED UNIT INTO TWO INDIVIDUAL BAGS AND HEAT SEAL ALL OPENINGS.
3. MINIMUM DOSE OF RADIATION NEEDED IS 25kGy AND GAMMA INDICATOR IS REQUIRED. COMPLETE ASSEMBLY AND DOUBLE BAGS TO BE STERILIZED.
4. THERE SHALL BE NO CHANGES IN THE ASSEMBLY PROCESS OR COMPOSITION OF MATERIALS WITHOUT PRIOR WRITTEN APPROVAL FROM MILLIPORE'S SUPPLIER QUALITY DEPARTMENT.
5. FOR ALL FIRST TIME BUILDS, MILLIPORE RESERVES THE RIGHT TO MODIFY PACKAGING WITHOUT NOTIFYING THE CUSTOMER.
6. WHEN MILLIPORE FILTERS AND/OR LYNK VALVES ARE INCLUDED IN THE ASSEMBLY, INCLUDE A COPY OF THE CERTIFICATE OF QUALITY FOR EACH WITH EACH SHIPMENT.
7. WHEN DUST COVERS ARE INCLUDED IN THE ASSEMBLY, SECURE WITH GAMMA COMPATIBLE TAPE, UNLESS OTHERWISE INDICATED, ALL CONNECTIONS SECURED WITH APPROPRIATE OETIKER CLAMP.

CUSTOMER APPROVAL:
SIGNED: _____
PRINT: _____
COMPANY: _____
DATE: _____

CUSTOMER APPROVAL: THIS DRAWING HAS BEEN PREPARED BASED ON INTERPRETATION OF CUSTOMER SPECIFICATIONS AND ASSUMES FULL RESPONSIBILITY THAT THE DESIGN CONFORMS WITH CUSTOMER REQUIREMENTS AND SIGNATURE INDICATES ACCEPTANCE OF RESPONSIBILITY AND APPROVAL OF THE DESIGN.

THIS DRAWING IS THE PROPERTY OF MILLIPORE CORPORATION AND IS SENT FOR CUSTOMER APPROVAL

MOBIUS by MILLIPORE

DATE: 12-MAY-11

DESIGNED BY: C. RAYE

PROFESSOR: NTS

DESIGNED BY: C. RAYE


DATE: 12-MAY-11

ADVANCING LIFE SCIENCE TOGETHER®
Research. Development. Production.

MCA0001L959A -CA 1

MOBIUS ASSEMBLY

Sterile Filtration




- Sterilizing filters for air filtration during integrity testing
- Flush bags and bags on vent/drain (operator safety)
- Gamma-irradiated single-use assembly (efficiency)
- Optimized hardware (ease of use)

Formulation and bulk preparation

Mixing, transfer and storage

Aseptic Alum Mixing Using Disposable Mixer



Final Formulation Filling Challenges

Limited Flexibility

- Varied product portfolio
- Fixed capacity
- Unable to respond to emergencies

Long Change-out Time

- Equipment set-up time
- CIP/SIP

Increasing Regulatory Qualification


- Product cross contamination
- Environmental control
- Operator safety

High capital equipment

Maintenance and spare parts costs


Operating Costs

- CIP (WFI)
- SIP
- Labor intensive
- Facilities




Vaccine Formulation & Filling

Highest level of product integrity and personnel protection




Changing the Paradigm

-Single-use Finish & Fill for vaccines



Design Considerations


Typical Final Filling Assembly: Good Design



Example: Final Formulation & Fill Finish Assembly




100L MIX Bag



Batch size: 100 L
Max. contact: 24 hours
Process temperature: RT
Product: Small molecule
Dosage: 1 mL syringe/daily


Sterile Transfer across a sterile barrier



Container/Bag approach

Lock by rotation (60°)

Open the double door



Case study: Single-Use Benefits

	Traditional	SU Solution
Clean and set-up	14 Hrs	<1 Hr
Cleaning validation	Extensive	Zero
Filling time	24 hrs	10 hrs
Average vials/hr	3,000	10,000
Aseptic connections	50	0
Operator Training	2 weeks	2 days
Equipment utilization	35%	82%
CAMPAIGN FILL TIME	36 Hrs	12 Hrs

Overview of single use products

Applications of single use technology in vaccine manufacturing


Validation of single use systems

Potential challenges with adoption of single use

Factors that may restrict use of disposables in biopharmaceutical manufacturing

Percent indicating "STRONGLY AGREE" or "Agree"

0% 10% 20% 30% 40% 50% 60% 70% 80%



8th Annual Report and Survey of Biopharmaceutical Manufacturing Capacity and Production – April 2011

Validation considerations

Incorporate QbD by selecting well qualified and safe materials (vendor selection)

- Components connections
- Integrity
- Sterilisation
- Packaging
- Shelf life
- Sterility

Vendor

Qualification data
Production controls

Certificate of Quality

Validation Guide

Supplier Audit

Manufacturer

Defined Product & Process Conditions



Performance

Risk assessment and qualification

- Chemical compatibility
- Extractable and leachable
- Impact on vaccine safety and efficacy
- Bioburden and endotoxin
- Stability studies

E&L PDA Definitions

Extractables


“Any chemical component that is removed from a material by the application of an artificial or exaggerated force (e.g., solvent, temperature or time).”

Determined under “worst-case” conditions following the Model Stream approach.

Leachables

“A chemical component that migrates from a contact surface into a drug product or process fluid during storage or normal use conditions.”

Determined with the product under normal processing/storage conditions.



Regulatory Agencies Expectations

*"Where there is **relevant risk**, the drug sponsor may have to **determine toxicity based on maximum dosage of potential leachables based on extractables data.**"*

*"**If there is no relevant risk** associated with the (material in question), **vendor data can be cross referenced** and a detailed justification for the applicability of these data and a justification for no additional testing should be submitted."*


Destry M. Sullivan - Senior Regulatory Review Officer, CBER
IBC's 7th International Single Use Applications for Biopharmaceutical Manufacturing
Conference, La Jolla, CA, June 14, 2010

Example: Final Formulation & Fill Finish Assembly



100L MIX
Bag

+



Batch size: 100 L
Max. contact: 24 hours
Process temperature: RT
Product: Small molecule
Dosage: 1 mL syringe/daily

Example: Collection of Vendor EXT Data

Component	TOC (mg C)
100L Mixer Bag	15.4
2L Bag	1.0
Connectors	0.36
Tubing	26.7
Tubing Manifold	47.1
Filter Capsule with Sterilizing Grade Membrane	30.1
TOTAL	120.7
Concentration of total extractables (mg/L = ug/ml)	3.01

**Bags Film
Vendor data**

**Connectors: Component
Vendor data**

**Tubing: Component
Vendor data**

Filter: Val. Guide

Assess Risk /
Criticality



Extractable
Evaluation



Assess Risk /
Patient
Safety



Leachables
or Additional
Testing



Assess Risk /
Patient
Safety

Points to Note on the Analysis / Approach

Did not require

- An in-house or consultant toxicologist
- The single-use system to be made and supplied
- The single-use system to be tested
- Specific analytical testing
- Parallels with container closure approach
- Anything other than a review of publically available documentation on extractables and leachables

HOWEVER it did rely on

- A qualified informed and experienced vendor
- An agreed final draft design
- An assigned person in the organization to be responsible
- A realistic timeline
- A multidisciplinary group in the organization

Risk Assessment Approach to identify Critical and Specific Service Needs



VENDOR

Process and Manufacturing
 Product and Patient Knowledge
 Internal Procedure and Controls
 Risk Tolerance
 Past Experience

Material/Component Knowledge
 Assembly Qualification and Design
 Manufacturing and Controls
 Assembly handling best practices
 Experience across many customer processes

Packaging Testing

Shelf Life

Sterilization Validation

Summary

- Use of single use technologies can quickly help increase operational flexibility and manufacturing capacity.
- Implementation of single use technologies is a multi-stage collaborative process between vendor and customer
- Great vendor support is critical to successful implementation and validation of SUS

Thank you !!!