




Tofflon

Presentation of QbD Application in Aseptic Filling Line

Jackson Zhao
2019-11-26

EXPERTISE IN PHARMACEUTICAL INDUSTRY



www.tofflon.com

© Shanghai Tofflon Science and Technology Co. Ltd. All Rights Reserved.

CONTENTS

1

Tofflon Filling Linkage Line in Injectable Solution System

2

QbD Application in Aseptic Filling Line

1 Filling Linkage Line in Injectable Solution System





Filling Linkage in Injectable System Solution

(1) Mini KUFILL

KUFILL—The Advanced Aseptic Processing System Integrated with Isolator

(2) KUFILL

(3) Flexible Manufacturing System Solution

(4) Vials Lyophilized System Solution

(5) Vials Liquid System Solution

(6) Ampoules Liquid System Solution

(7) Vials Powder Filling System Solution



Mini KUFILL




KUFILL




Flexible Manufacturing System Solution




Vials Lyophilized System Solution



Vials Liquid System Solution



Ampoules Liquid System Solution



Vials Powder Filling System Solution





2 Quality by Design

Quality by Design

Product Quality Lifecycle Implementation (PQLI)




Quality by Design: A systematic approach to development that begins with predefined objectives and emphasizes **product and process understanding and process control**, based on sound science and quality risk management.



Quality by Design

A Systematic Approach



- Ishikawa • Target the product profile
- Capability • Determine Critical Quality Attributes (CQAs)
- FMEA • Link input material attributes and process parameters to CQAs and perform risk assessment
- Pareto • Develop a design space
- DOE • Design and implement a control strategy
- Manage product lifecycle, including continual improvement

Design Management

Aseptic Filling Line




Design Team 01

Experienced engineers and expert advisors

02 Design Tools

Pro/E and Inventor/CAD used for 2D/3D design, PDM and OA used for drawings management.



Design Plan 04


A precise plan, plan management and quick response for each project according to the design flow chart.

03 Design Reference

GMP/GEP/GAMP5 guidelines and 21 CFR Part 211 and Part 11, etc.


Design Management

Design Reference




Design Management

Design Plan




Design Management

Lifecycle



Design Management

Review/Approve/Release




Design Management

Technical Clarification



Design Management


Change Control



- Here we introduce the **QbD** concept to the **Quality Control of the Aseptic Filling Line**
- **QbD Three main tools:**
 - ❖ **Design Of Experiments(DoE)**
 - ❖ **Failure Mode and Effects Analysis(FMEA)**
 - ❖ **Process Analytical Technology(PAT)**

The First QbD tool: DoE

QbD-DoE



To build a design space and calculate parameter of the whole filling line accurately and ensure stability, reliability and safety of equipment

Product Performance (FS DS HDS SDS)

Process Performance (IQ OQ PQ Test Protocol)

Good SOPs and perform SOPs during fabricating process

Product Design

Technical Data



Technical Data For Washing Machine

Tofflon®	No. : TFL/	
	Ver. : B	Page : 2 / 11
Title : TD for Washing machine		
Owner : Shanghai Tofflon Science & Technology Co., Ltd.		

INDEX

1. Processing range and performance.....	3
1.1. Main parameter.....	3
1.2. Vials output.....	3
1.3. Dimensions and weight.....	4
1.4. Utility Requirements.....	4
2. Calculation.....	4
2.1. Parameter analysis and confirmation.....	4
2.1.1. The operating cycle of the equipment.....	4
2.1.2. The number of needles.....	5
2.1.3. The size of pitch.....	5
2.2. Calculation of the transmission ratio of washing machine.....	5
2.2.1. Calculation of the needles action time.....	5
2.2.2. Calculation the spraying time of needle.....	5
2.2.3. Selection of number of spraying needles.....	5
2.3. Transmission design.....	6
2.3.1. Drive mechanism.....	6
2.3.2. Transmission chain.....	6
2.4. Pitch design.....	7
2.5. Cam design.....	8
3. Utility supply.....	9
3.1. The consumption of WFI and CA.....	9
3.2. Thermal radiation analysis.....	9
4. Washing machine Process description.....	10
4.1. main operation modes.....	10
4.2. parameters management.....	11
4.3. recipe management.....	12
4.4. Drain model.....	13
5. Revision record.....	13

Technical Data For Sterilizing Depyrogenation Tunnel

Tofflon®	No. : TFL/	
	Ver. : C	Page : 2 / 15
Title : TD for Tunnel		
Owner : Shanghai Tofflon Science & Technology Co., Ltd.		

1. Output table.....	3
1.1. Vials output.....	3
1.2. Dimensions and weight.....	4
1.3. Utility requirements.....	4
2. Calculation for length of heating zone.....	4
2.1. Laws and regulations.....	4
2.2. Calculation for dwell time.....	4
2.3. Calculation for speed for belt.....	6
3. Calculation for the power of the heater.....	7
3.1. Regular empirical method.....	7
3.2. Theoretical method for heater calculation.....	7
4. Select the fan and calculate air volume.....	9
4.1. Air volume for heating zone.....	9
4.2. Air pressure for heating zone.....	9
5. Calculation of cooling section.....	9
5.1. Heat carried away by bottle cooling-E1.....	9
5.2. Heat carried away by belt cooling -E2.....	10
5.3. Select the number of filter in the cooling zone.....	10
6. The heat-exchange for chilly water cooling.....	11
6.1. Operating conditions:.....	11
6.2. Based on above parameter conditions,.....	11
6.3. Chilled water flow(L).....	11
6.4. Cooler ventilation area(SC).....	11
6.5. Average heat transfer temperature difference (Δt_m).....	11
6.6. Heat exchanger area(F).....	11
6.7. Select the number of rows in the heat exchanger.....	11
7. Tunnel Process description.....	12
7.1. HEAP pre-heating process.....	12
7.2. Main operation modes.....	12
7.3. Sterilization of cooling zone.....	14
7.4. Air pressure balance control.....	15
8. Revision record.....	17

Technical Data For Filling and Stoppering Machine

Tofflon®	No. : TFL/	
	Ver. : G	Page : 2 / 23
Title : TD for Filling and Stoppering Machine		
Owner : Shanghai Tofflon Science & Technology Co., Ltd.		

INDEX

1. Processing range and Performance.....	3
1.1. Processing range.....	3
1.2. Performance.....	3
2. Calculation.....	3
2.1. Parameter analysis and confirmation.....	3
2.1.1. Weighing system analysis.....	4
2.1.2. Vials conveying speed analysis.....	5
2.1.3. Filling process analysis.....	5
2.1.4. Summary.....	5
2.2. Timing and phase design.....	6
2.3. Transmission design.....	6
2.3.1. Drive chain.....	7
2.3.2. Confirmation of bottle feeding capacity.....	9
2.3.3. Feeding system.....	11
2.3.4. Filling related structure.....	13
2.3.4.5. Cam design.....	17
2.4. Filling process.....	18
2.4.1. Filling process.....	18
2.4.2. Pump selection.....	19
2.5. The process of IPC.....	19
2.5.1. The process of IPC.....	19
2.5.2. Load cell selection.....	20
3. Production capacity.....	20
4. Design output and conclusion.....	22
5. Revision record.....	23

Product Design

FS/DS/HDS/SDS/DQ



Tofflon

FUNCTIONAL DESIGN SPECIFICATION 功能设计规格书

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	NO./编号: 2016-028LF
	MODEL/型号: FFVLP50/01-01
FUNCTIONAL DESIGN SPECIFICATION FOR FILLING AND STOPPERING MACHINE 灌装压塞机功能设计规格书	

Tofflon

HARDWARE DESIGN SPECIFICATION 硬件设计规格书

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	NO./编号: 2016-028LF
	MODEL/型号: FFVLP50/01-01
HARDWARE DESIGN SPECIFICATION FOR FILLING AND STOPPERING MACHINE 灌装压塞机硬件设计规格书	

Tofflon

SOFTWARE DESIGN SPECIFICATION 软件设计规格书

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	NO./编号: 2016-028LF
	MODEL/型号: FFVLP50/01-01
SOFTWARE DESIGN SPECIFICATION FOR FILLING AND STOPPERING MACHINE 灌装压塞机软件设计规格书	

Tofflon

DESIGN QUALIFICATION 设计确认

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	NO./编号: 2016-028LF
	MODEL/型号: FFVLP50/01-01
DESIGN QUALIFICATION FOR FILLING AND STOPPERING MACHINE 灌装压塞机设计确认	

Process Design-Manufacturing Quality Control

Manufacturing Quality Control



Good SOPs



and perform SOPs during fabricating process

IQC (Incoming Quality Control)
IPQC (InPut Process Quality Control)
FQC (Finish or Final Quality Control)
OQC (Outgoing Quality Control)




The Second QbD tool: FMEA

ICH Q9 Quality Risk Management



Quality by Design

FMEA



S:Severity O:Occurrence D: Detection

Quality by Design

SIA/CCA

Tofflon



Tofflon

SYSTEM IMPACT ASSESSMENT

系统影响性评估

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	PROJECT NO./项目编号: TFL-IPS-2016-034
	CONTRACT NO./合同编号: ZMCBE1-1609005
SYSTEM IMPACT ASSESSMENT 系统影响性评估	

COMPONENT CRITICALITY ASSESSMENT

部件关键性评估

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	NO./编号: 2016-028LF
	MODEL/型号: FFVLP50/01-01
COMPONENT CRITICALITY ASSESSMENT FOR FILLING AND STOPPERING MACHINE 灌装压塞机部件关键性评估	

RISK ASSESSMENT

风险评估

SUPPLIER NAME/供应商: Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称: Zhejiang Medicine Co., Ltd Xinchang Pharmaceutical Factory 浙江医药股份有限公司新昌制药厂
	NO./编号: 2016-028LF
	MODEL/型号: FFVLP50/01-01
RISK ASSESSMENT FOR FILLING AND STOPPERING MACHINE 灌装压塞机风险评估	

Quality by Design

RA



Tofflon	RA	Doc. No. 文件编号	2016-028LF-RA
		Model 型号	FFVLP50/01-01
Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	Risk Assessment for Filling and Stoppering Machine 灌装压塞机风险评估	Version 版本号	V 1.0
		Page 页码	14 of 40

9. Risk Assessment Implementation and Record 风险分析执行及记录

Number 编号	Component/ Function 部件/功能	Potential failure mode 潜在失效模式	Potential failure consequences 潜在失效后果	S	Potential failure cause 潜在失效原因	O	D	R P N	Control Measure 控制措施	Qualificat ion Activity 确认活动	P S	P O	P D	PR PN	Whethe r accept 是否接受	Execution confirmati on 执行情况 确认
YS251	Minimum load sensor 缺瓶传感器	Error display 显示错误	Affect the normal operation of machine 影响设备正常运行	5	Program failure or sensor failure 程序故障或传感器故障	1	1	5	Alarm confirmation. 进行报警测试	1.OQ	5	1	1	5	Y 是	
YS006	Minimum amount of warning sensor 最小量预警 传感器	Error display 显示错误	Affect the normal operation of machine 影响设备正常运行	5	Program failure or sensor failure 程序故障或传感器故障	1	1	5	Alarm confirmation. 进行报警测试	1.OQ	5	1	1	5	Y 是	
YS005	Minimum load sensor 缺瓶/停机传 感器	Error display 显示错误	Affect the normal operation of machine 影响设备正常运行	5	Program failure or sensor failure 程序故障或传感器故障	1	1	5	Alarm confirmation. 进行报警测试	1..OQ	5	1	1	5	Y 是	
M178	Main motor 主电机	Transports insufficient capacity	Affect the normal operation of machine, eg,	5	Improper motor selection 电机选型不当	1	3	15	Correct installation confirmation 正确安装确认	1.IQ	5	1	1	5	Y 是	

The Third QbD tool: PAT

- Here We Use the Following List of Instruments to Ensure the Quality of Aseptic Filling Line.

Manufacturing Quality Control

IQC



S.N.	IQC SOPs and Check Lists
01	IQC SOP
	SOP
	Flow chart
02	Key parts inspection SOP
	Weighing module
	Manifold
	Filling needles
	buffer tank
	Stoppering air chamber
	Rejection air chamber

03	Packing materials management
	SOP
	Flow chart
	Check list
04	Materials management
	Stainless steel management SOP
	Mental parts (Stainless steel) and non-mental parts SOP
05	Unqualified parts handling
06	Supplier management
	Flow chart and requirement for inspection at supplier site
	Inspection by supplier SOP
	Quality assurance agreement
	Supplier processing requirement

Manufacturing Quality Control

IQC-Material Certification Management-URS Review



N°	Requirement description	Criticality (Yes/No)
Generalities		
U.5.2.1	Electrical panels must be made of steel.	No ⁽²⁾
U.5.2.2	All metal parts in direct product contact or in contact with the stoppers must be made of stainless steel 316L (1.4404). Material Certificate of Conformity 3.1 must be provided. The surface finish / rugosity must be lower or equal to 0.6µm.	Yes
U.5.2.3	All non-metal parts must be in conformity with the USP class VI (certificates must be provided)	Yes
U.5.2.4	QSMP sheets must be provided for all parts and elements (hoses) in direct product contact.	Yes
U.5.2.5	All metal parts in direct contact with process compressed air must be in stainless steel 316L (1.4404). Material Certificate of Conformity 3.1 must be provided. The surface finish / rugosity must be lower or equal to 0.8 µm.	Yes

Manufacturing Quality Control

IQC-Material Certification Management



□ Material Management Flow Chart



Manufacturing Quality Control

IQC-Material Certification Management



Material introduction inside the RABS Cabinet (Conveyor and Filling machine)

No.	Item/Name	MOC/Spec./Desc.	Desc. For Contact	Location	Reference(certification)
1	Infeed starwheel	POM	Contact with external of the vial	Infeed system	FDA / ISO
2	Filling needles	SS 316L	contact with Product	Filling system	GB/T 20878-2007
3	Manifolds and pipes	SS 316L	contact with Product	Filling system	ASME BPE
4	Filling hoses	Silicone	contact with Product	Filling system	FDA / USP Class VI
5	Filling needles for gas flushing	SS 316L	contact with Product	Filling system	GB/T 20878-2007
6	Manifolds and pipes for gas flushing	SS 316L	contact with Product	Filling system	ASME BPE
7	Filling hoses for gas flushing	Silicone	contact with Product	Filling system	FDA / USP Class VI
8	Body of diaphragm valve	SS 316L	contact with Product	Filling system	EN10204 3.1
9	Diaphragm of diaphragm valve	EPDM/PTFE	contact with Product	Filling system	FDA/BSE/TSE/USP Class VI
10	Buffer tank	SS 316L	contact with stopper	Filling system	GB24511-2009 / ASTM-A240-316L / GB/T 25198-2010
11	Pipelines for WFI, sterile compressed air and pure steam	SS 316L	contact with Product	CIP/SIP system	ASME BPE
12	Pipelines for Product	SS 316L	contact with Product	CIP/SIP system	ASME BPE
13	Filter housing	SS 316L	contact with Product	CIP/SIP system	Statement / EN10204 3.1
14	Body of diaphragm valve	SS 316L	contact with Product	CIP/SIP system	EN10204 3.1
15	Diaphragm of diaphragm valve	EPDM/PTFE	contact with Product	CIP/SIP system	FDA/BSE/TSE/USP Class VI
16	Sealing ring	PTFE	contact with Product	CIP/SIP system	FDA / USP Class VI
17	Bowl and chute	SS 316L	contact with stopper	Stoppering system	ASTM A240 ASME SA-240 / EN10204 3.1
18	Stoppering plate	SS 316L	contact with stopper	Stoppering system	ASTM A-240M/480M / EN10204 3.1
19	Accelerative wheel	POM/SS304	Contact with external of the vial	Main transmission system	FDA / ISO
20	Transfer toothed plate	SS 304	Contact with external of the vial	Main transmission system	EN10028-7:2007 / EN10204 3.1
21	Rejection starwheel	POM	Contact with external of the vial	Rejection system	FDA / ISO
22	Discharging starwheel	POM	Contact with external of the vial	Outfeed system	FDA / ISO
23	Discharge belt	POM	Contact with external of the vial	Outfeed system	FDA / ISO
24	Guide	POM /PE	Contact with external of the vial	Outfeed system	FDA / ISO
25					

Manufacturing Quality Control

IQC-Material Certification Management



Tofflon	DQ	Doc. No.	2017-017LFe-DQ	
		Model	FFVLP40/01-01	
		Version	V 1.0	
		Page	61 of 73	
Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	Design Qualification for Filling and Stoppering Machine 灌装压塞机设计确认	Design Qualification for Filling and Stoppering Machine 灌装压塞机设计确认	Design 设计情况	Yes/No 是否符合要求
10.6. Materials and Surface Finishes Qualification 材料和表面 <u>Purpose 目的</u> Verify whether materials of surface that direct or indirect cGMP/URS/TS requirements. 确认与产品直接或间接接触的材料是否符合 cGMP/URS/TS 要求。 <u>Acceptance Criteria 可接受标准</u> All materials of surface that direct or indirect cGMP/URS/TS requirements. 所有与产品直接或间接接触的材料均应符合 cGMP/URS/TS 要求。		Design	Yes (✓) No ()	
		Design	Yes (✓) No ()	
		Design	Yes (✓) No ()	
		Design	Yes (✓) No ()	
TS4.4.12 TS4.4.13	Starwheel 星轮	POM	POM	Yes (✓) No ()
TS4.4.11	V-block 输入块	POM	POM	Yes (✓) No ()
TS4.4.14	Belt 网带	POM	POM	Yes (✓) No ()

EXPERTISE IN PHARMACEUTICAL INDUSTRY
专业技术服务于制药工业

Yuki. Wu 2018.10.17

Tofflon	DQ	Doc. No.	2017-017LFe-DQ	
		Model	FFVLP40/01-01	
		Version	V 1.0	
		Page	63 of 73	
Shanghai Tofflon Science & Technology Co., Ltd. 上海东富龙科技股份有限公司	Design Qualification for Filling and Stoppering Machine 灌装压塞机设计确认	Design Qualification for Filling and Stoppering Machine 灌装压塞机设计确认	Design 设计情况	Yes/No 是否符合要求
Conclusion/Remarks 结论备注 Remark 备注: The acceptance criteria is fulfilled. no deviation is found. 印+批 2018.10.17 Conclusion 结论 Conformity 符合性 Yes/是 (✓) No/否 () Deviation No. 偏差编号 N/A Executed by / Date: 执行人/日期 Yuki. Wu 2018.10.17 印+批 2018.10.17 Reviewed by / Date: 审核人/日期 Anne Wilson 2018.10.17 Oleg Finkov 2018.10.17				

EXPERTISE IN PHARMACEUTICAL INDUSTRY
专业技术服务于制药工业

www.tofflon.com

Manufacturing Quality Control

IQC-Material Certification Management



Tofflon

MATERIAL OF CERTIFICATE

材质证明

SUPPLIER NAME/供应商: Shanghai Tofflon Science&Technology Co.,Ltd. 上海东富龙科技股份有限公司	CUSTOMER NAME/客户名称 *****
	NO./编号: 2017-017LFe
	MODEL/型号: FFVLP40/01-01
MATERIAL OF CERTIFICATE FOR FILLING AND STOPPERING MACHINE 灌装压塞机材质证明	

EXPERTISE IN PHARMACEUTICAL INDUSTRY

ADD: 139, DUHUI RD., SHANGHAI -201108 CHINA
TEL: +86 21 64909996
FAX: +86 21 64908881
Http: // www.tofflon.com



Tofflon	MOC	Doc. No. 文件编号	2017-017LFe-MOC
		Model 型号	FFVLP40/01-01
		Version 版本号	V 1.0
Shanghai Tofflon Science&Technology Co.,Ltd. 上海东富龙科技股份有限公司	Material Certificate for Filling and Shopping Machine 灌瓶压盖机材质证明	Page	1 of 3

Material of Certificate List 材质证明清单

No. 序号	Component 部件	Material 材质	Critical or Non- critical (C or NC) 关键或非关键	Remark 备注
1	Filling needles 灌装针	SS316L	C	
2	Manifold 分配管	SS316L	C	
3	Pipeline 管路	SS316L	C	
4	Filling tubes 灌装软管	Silicone tube 硅胶管	C	
5	Filling pump 灌装泵	SS 316L	C	
6	Buffer tank 缓冲罐	SS316L	C	
7	Body of diaphragm valve 隔膜阀本体	SS 316L	C	
8	Diaphragm of diaphragm valve 隔膜阀膜片	EPDM	C	
9	The stopper feed bowl and feed tracks 胶塞料斗和供料轨道	SS 316 L	C	
10	Turn table plate 转盘板	SS 304	C	

Manufacturing Quality Control

IQC-Material Certification Management



MATERIAL CERTIFICATE

Customer: 东富龙 Certificate No.: 180622-412
Customer Order No.: 3215TFLF180517001 Date of Certificate: 2018-05-22
Part No.: FFF17019201400901 Production Batch No.: N/A
Part Description: 10 头分配管 Dosing vessel.

Raw Material Specification

Heat Number	Material Specification	Material suppliers	Material Description	C	Si	Mn	P	S	Ni	Cr	N	Mo
160325W07	Ø55	浦青	316L	0.020	0.32	0.96	0.037	0.0066	10.08	16.24		2.02
17ZT0107	Ø12	浙泰不锈钢	316L	0.016	0.30	0.92	0.036	0.008	10.08	16.37		2.01
51076	Ø38.1±1.65	凌士通	316L	0.020	0.48	0.88	0.035	0.007	10.00	16.60	0.042	2.04

Mechanical Test

Heat Number	Yield Rp 0.2		Yield Rp 1.0		Tensile Strength		Hardness	Reduction	Elongation	Impact Test
	N/mm2	PSI	N/mm2	PSI	N/mm2	PSI	HRB/HB	%	%	20°C-J
160325W07	-	-	-	-	-	-	-	-	-	-
17ZT0107	-	-	-	-	-	-	-	-	-	-
51076	295	-	-	-	625	-	82	-	56	-

Mechanical Test

Heat Number	Eddy Current Test	Visual & Dimensional Test	Flaring Test	Flattening Test	Intergranular Corrosion Test	Material Identification Test
160325W07	OK	OK	OK	OK	-	OK
17ZT0107	OK	OK	OK	OK	-	OK
51076	OK	OK	OK	OK	-	OK

We certify that this information is a true representation of the data that has been furnished by our raw material suppliers. We have no knowledge of any mercury of low melting contamination.

Examined according to quality system ISO 9001:2008 by QAI. Certificate number: QAIC / CN / 121530-A

This certified Mill Test Report is computer generated and meets EN 10204 3.1

This Certificate was made by use of a computer system and is valid without signature.



ASME BPE Compliant



Expert Precision Material Engineering & Technical Ltd.

Manufacturing Quality Control

IQC-Material Certification Management



Manufacturer:
34 rue du Moulin aux Aulnaies
89120 Charny-Orée-de-Puisaye, France

Certificate of Conformance

Part Number/ Revision:	760268	Customer Part Number/ Revision:	84-311-032 Cdc rev.02
Description:	T. VERSILIC 3.2X6.4- 25M MARQUE - DOUBLE SACS		
Lot Number:	1959481	Lot Quantity:	4000 m
Date of Manufacture: (DD/MMM/YYYY)	24/04/2018	Expiration Date: (DD/MMM/YYYY)	24/10/2020
Post processing Run Number: (Refer to the attached Certificate of Processing for Additional Detail)		N/A	

Filling Tubes:
Silicone

We certify the material listed above conforms in full with the following specifications:

All items have been manufactured, inspected, tested, and accepted in accordance with our Quality Management System certification.

All materials and notes indicated on the purchase order, drawing, specifications, quality assurance requirements, Versilic® Regulatory Information Overview (RIO), Revision 0, or other applicable approved documents effective on the date of manufacture. The ASTM plaques of Versilic® silicone meet the requirements of the USP, Biological Test for Plastics, Class VI.

All materials and processes used in manufacturing conform to the materials and/or manufacturing specifications and notes indicated on the purchase order, drawing, specifications, quality assurance requirements, Versilic® Regulatory Information Overview (RIO), Revision 0, or other applicable approved documents effective on the date of manufacture. The ASTM plaques of Versilic® silicone meet the requirements of the USP, Biological Test for Plastics, Class VI.

Gobain's sole liability shall be to replace any product not in conformance with the specification and requirements of the buyer.



Quality Approval:	Christelle PEDROSO	Date:	27/04/2018
-------------------	--------------------	-------	------------

Manufacturing Quality Control

IPQC-Mechanical Installation






TFL/RD7.5-01-22A-LF.

Tofflon		灌装机装配检查表			
客户名称:		机器型号:			
机器编号:		WB S(组件):		灌装机关键部位:	
检查项目	参考值	实际值	执行人	过板检	
 <p>测量图 1 机器各段的长度, 对角线, 保证符合设计尺寸, 误差 $\leq 1\text{MM}$。</p>	$\leq 1\text{MM}$				
 <p>测量图 2 机器主体水平, 应每米 $\leq 0.5\text{MM}$。〈装线包完成之前〉。</p>	$\leq 0.5\text{MM}$				
 <p>图 3 所示, 测量各类型机侧边错位, 应 $\leq 0.2\text{MM}$。</p>	$\leq 0.2\text{MM}$				
 <p>图 4 所示, 测量各类型机顶部错位, 应 $\leq 0.2\text{MM}$。</p>	$\leq 0.2\text{MM}$				

1/9

TFL/RD7.5-01-22A-LF.

Tofflon		灌装机装配检查表			
客户名称:		机器型号:			
机器编号:		WB S(组件):		灌装机关键部位:	
检查项目	参考值	实际值	执行人	过板检	
 <p>转盘安装完后检查其转动, 其转动 $\leq 0.5\text{MM}$。轴向跳动 $\leq 0.15\text{MM}$。</p>	$\leq 0.5\text{MM}$				
 <p>转盘与排液阀用隙均匀一致, 应 $\leq 2\text{MM}$。</p>	$\leq 2\text{MM}$				
 <p>排液阀与转盘 $0.1-0.2\text{MM}$。</p>	$0.1-0.2\text{MM}$				

2/9

TFL/RD7.5-01-22A-LF.

Tofflon		灌装机装配检查表			
客户名称:		机器型号:			
机器编号:		WB S(组件):		灌装机关键部位:	
检查项目	参考值	实际值	执行人	过板检	
 <p>小转盘与进料量轮等轴同轴度 $\leq 1.5\text{MM}$。不可偏转转盘。 2. 转盘与排液阀 $0.1-0.2\text{MM}$。</p>	$\leq 1.5\text{MM}$				N/A
 <p>如图, 转盘量轮与排液量轮转动良好, 互相之间无干涉现象, 各量轮轴与轴跳动与轴肉跳动, 应 $\leq 0.05\text{MM}$。</p>	$\leq 0.05\text{MM}$				

3/9




100





Toflon

电器柜装配自检记录表。

客户名称:		机型:			
项目编号:					

检查项目:	检查项:	标准:	执行人:	复核人:
	检查Y20侧星型电源线是否有与零地短路。大地可算零按《独立电零》。			
	接线保证零地线。一个桩脚不允许有三个以上的零地线。			
	检查右线是否有强电分开。接强电在左。弱电在右分开。			

Tofflon		电器柜装配自检记录表			
客户名称		机型			
项目编号					
检查项目		检查项	结果	执行人	复查人
		检查 DIN导轨 双轴连接块 对应接口 标签与角 量电机组 上标签一 致 Y/N			
		复核人员 检验确认 后,在图纸 上用荧光 笔做标记 Y/N			
		配线架接 线检查后 应,安装电 源柜,确认 电箱柜与 底座接文 架是否牢 固 Y/N			

Tofflon		电器柜装配自检记录表		
客户名称		机型		
项目编号				
检查项目 		检查项目 检查欧美项目是否有出现国际插座(外插电源) Y/N		
		检查项目 欧美项目色标开关必须需外亮 Y/N		
		检查欧美项目 三色灯是否有指示标记 Y/N		

Manufacturing Quality Control

FQC-Test Protocols



S.N.	FQC Check Lists of Washing Machine
01	Preparation before starting up
02	Alarm Test
03	IO Test
04	P&ID Drawing Verification
05	Qualification of Safety Equipment
06	Check HMI Screen
07	Motors Verification
08	Welding Verification
09	Check of Control System
10	Layout Arrangement Verification
11	Qualification of Washing Process
12	System Access
13	System Operational Test
14	Noise Permitted Test

S.N.	FQC Check Lists of Tunnel
01	Preparation before starting up
02	Alarm Test
03	IO Test
04	P&ID Drawing Verification
05	Qualification of Safety Equipment
06	Check HMI Screen
07	Fan Verification
08	Check of Control System
09	Layout Arrangement Verification
10	System Access
11	Belt Speed Test
12	Empty Chamber Temperature Distribution Test
13	Full Load Temperature Penetration Test
14	Surface Temperature Test
15	System Operational Test
16	Noise Permitted Test

S.N.	FQC Check Lists of Filling Machine
01	Preparation before starting up
02	Alarm Test
03	IO Test
04	P&ID Drawing Verification
05	Qualification of Safety Equipment
06	Check HMI Screen
07	Motors Verification
08	Welding Verification
09	Check of Control System
10	Layout Arrangement Verification
11	System Access
12	Filling Accuracy Test
13	System Operational Test
14	Counting Test
15	Noise Permitted Test

Manufacturing Quality Control

FQC-Test Protocols



Tofflon	Filling Accuracy Test 灌装精度测试	Doc. No. 文件编号	SOP-053
		Version 版本号	V1.0
Shanghai Tofflon Science & Technology Co., Ltd. 上海东龙科技股份有限公司	Filling and Stoppering Machine 灌装压盖机	Page 页码	3 of 7
		Issued Date 发布日期	2019.04.30

Tofflon	Filling Accuracy Test 灌装精度测试	Doc No 文件编号	SOP-053
		Version 版本号	V1.0
		Page 页码	4 of 7
Shanghai Tofflon Science & Technology Co., Ltd. 上海东龙科技股份有限公司	Filling and Stoppering Machine 灌装压盖机	Issued Date 发布日期	2019.04.30

Tofflon	Alarm Test 报警测试	Doc No 文件编号	
		Model 型号	
		Version 版本号	V 3.4
		Page 页码	3
Shanghai Tofflon Science&Technology Co Ltd 上海东龙科技股份有限公司	Filling and Stoppering Machine 灌装压盖机		

Filling Accuracy Test
灌装精度测试

Test result / 测试结果			
The filling pump filling accuracy is in the range of $\pm 1\%$. 灌装泵灌装精度在 $\pm 1\%$ 之内。			Yes () / No ()
The surrounding test conditions / 测试条件			
Starting time 测试开始时间	T1	Ending time 测试结束时间	T1
	T2		T2
	T3		T3
Containers size 容器规格	Room temperature (°C) / humidity (%) 房间温湿度		
Filling system 灌装系统	<input type="checkbox"/> Piston pump 柱塞泵灌装系统 <input type="checkbox"/> Peristaltic pump 蠕动泵灌装系统 <input type="checkbox"/> Auger dosing 螺杆计量(粉末) <input type="checkbox"/> Others 其他		
Filling product 灌装产品	Volume of piston pumps 柱塞泵体积		
Density 密度	Ø tube inside of peristaltic pump 蠕动泵的管内径 Ø		
Production speed 生产速度	Filling needle inside (Ø) 灌装针内径(Ø)		
Ids of pumps (Only for piston pump) / 泵编号 (柱塞泵适用)			
Pump 1/泵 1		Pump 2/泵 2	
Pump 3/泵 3		Pump 4/泵 4	
Pump 5/泵 5		Pump 6/泵 6	
Pump 7/泵 7		Pump 8/泵 8	
Pump 9/泵 9		Pump 10/泵 10	
Pump 11/泵 11		Pump 12/泵 12	

Test 1 测试 1

No. 泵号	Pump 1 1号泵	Pump 2 2号泵	Pump 3 3号泵	Pump 4 4号泵	Pump 5 5号泵	Pump 6 6号泵	Pump 7 7号泵	Pump 8 8号泵	Pump 9 9号泵	Pump 10 10号泵	Pump 11 11号泵	Pump 12 12号泵
Weight 重量												
Time 取料时间												
Maximum Value 最大值												
Minimum Value 最小值												
Yes/No 是否符合 要求	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()	Yes () No ()
Acceptance Value 接受标准	Filling Set Value 设定容量: <u> </u> ml			Filling Accuracy 灌装精度: <u> </u> %			Filling Water Acceptability Limit 卸料量合格范围: <u> </u> ml					

If necessary, this page can be copied. 如有需要, 此页可复制

Alarm Test
报警测试

Test Description 测试描述	Test Procedure 测试方式	Expected result 预期结果	Yes/No 是否符合要求
Emergency stop alarm 急停报警	Press the relevant emergency stop button 按下相应急停按钮	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止运行, 显示报警信息	Yes () No ()
Vacuum low pressure 真空低报警	Disconnect vacuum 断开进口真空	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止运行, 显示报警信息	Yes () No ()
CA low Pressure 压缩空气低报警 (如有)	Disconnect CA 断开进口压缩空气	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止运行, 显示报警信息	Yes () No ()
Vials infed belt minimum load 进瓶侧缺瓶报警	During production, prevent vials into the machine manually 在生产过程中, 人为的阻止瓶子进入	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止, 显示报警信息	Yes () No ()
Outfeed vials blocking alarm 出瓶堵塞报警	During production, manually vials are filled on outfeed belt 在生产过程中, 人为将出瓶通道堵塞瓶子	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止, 显示报警信息	Yes () No ()
Reject system is full 剔除满瓶报警	During production, manually block the sensor 在生产过程中, 用手堵住相应传感器	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止, 显示报警信息	Yes () No ()
Stopper feed channel minimum load 胶塞通道缺塞报警	During production, manually prevent stoppers into the feed channel 在生产过程中, 人为阻止胶塞进入胶塞通道	The lamp is on, the horn beeps, and the machine stop, display the alarm information. 报警灯亮, 蜂鸣器响, 机器停止, 显示报警信息	Yes () No ()

Manufacturing Quality Control

Test Instruments List



- Alloy analyzer: Innov system, U.S. brand
- Ultrasonic thickness inspector: imported from Japan
- Roughness inspector: imported from Japan
- Three-coordinate Detector: imported from France
- Video Measuring System: TaiShuo
- Temperature validation system: imported from UK
- Temperature calibrator: imported from UK
- Pressure calibrator: imported from UK
- Temperature validation system: GE KAYE, U.S. brand
- Temperature calibrator: GE KAYE, U.S. brand
- Signal generator: VICTOR
- Standard resistance box: ZX25a type
- Electron microscope: UK brand
- Endoscope detector: Olympus



Tofflon Vision: Smart Pharma Factory Builder

**Warmly Welcome
To Tofflon**

