

## Poll 1

Talk to our experts  
about the mRNA  
related processes  
at

**Raman.Rani@  
Sartorius.com**

## Poll 2

Talk to our experts  
about the mRNA  
related processes  
at

**Raman.Rani@  
Sartorius.com**

# SARTORIUS

Simplifying Progress

**Accelerating mRNA process using  
scalable, flexible and reliable  
solutions**

DCVMN, 12 April 2023




# Agenda

- Introduction
- Upstream – *in vitro* transcription
- Downstream Solutions
- Formulation Development & Storage
- Data Analytics



# Understanding mRNA Structure and Production



- Single-stranded RNA
- 5' Cap and poly-A tail are required for successful protein expression in cells.
- *In Vitro* transcription produces RNA from a DNA template, often a plasmid DNA
- 5' Cap can be added co-transcriptionally (during IVT), or post-transcriptionally
- Poly-A tail can be encoded in the DNA template, or added enzymatically after IVT

Good understanding of IVT reactions is fundamental to maximise productivity and document the purity.

# The Benefits of mRNA in Vaccines

Use the patient's cells as a manufacturing facility and patient creates its own vaccine.

mRNA  
technology...

01

**Reduce timeframes:** development and manufacturing (including variants update)

02

**Safety:** Cell free and virus free production.

03

**Footprint:** no seed train, smaller reactors...

mRNA : messenger Ribonucleic Acids

# The Challenges of mRNA Technology

mRNA  
technology...

01

**Doses** is dependent of application ( infectious diseases vs cancer)

02

**Size, sequence and type** of construct.(i.e new constructs such as saRNA are very long).

03


**Cost per dose** due to raw material (enzymes, capping...)

03

**New field, low expertise, no platform.**


mRNA : messenger Ribonucleic Acids

# mRNA: doses challenge...





# (m)RNA: various constructs to overcome some challenges BUT is challenging the process...



## messenger RNA (mRNA)


- Based on 'standard' cell mRNA

### Pros

- ✓ Short transcript (easier to manufacture)
- ✓ Lots of human trial data

### Cons

- ✗ Non-replicating
- ✗ One Gol per transcript



## self-amplifying RNA (saRNA)


- Based on family of self-replicating alphavirus (Built-in RNA polymerase (viral replicon))

### Pros

- ✓ Requires **lower dosage** for equivalent efficacy (non-vaccine applications)
- ✓ More durable expression than non-replicating mRNA

### Cons

- ✗ Longer than traditional mRNA (difficult to manufacture)
- ✗ Viral origin - greater immunogenicity



## trans-amplifying RNA (taRNA)


- Like saRNA, but split (trans) replicon/Gol

### Pros

- ✓ Requires **lower dosage** for equivalent efficacy (non-vaccine applications)
- ✓ Shorter than saRNA (easier to manufacture)

### Cons

- ✗ Requires delivery of two separate transcripts (replicon & Gol)
- ✗ Viral origin - greater immunogenicity



## circular RNA (circRNA)

- synthetic form of RNA

### Pros

- ✓ Durable expression & lower immunogenicity (closed configuration, non-vaccine applications)
- ✓ Lack of cap & polyA tail (easier manufacture)

### Cons

- ✗ Lack of clinical validation (still in pre-clinical stage)
- ✗ Manufacturing scalability not established (esp for longer RNA sequences)

What does it mean for the process/facility?

# Key Drivers for mRNA Manufacturing

## Current Technology

- Use plasmid for DNA template
- In vitro transcription
  - Enzymatic capping
  - Co-transcriptional capping
- LNP-Based Formulation
- Downstream needs to be adjusted to fit different capping methods

## Pain Point

- High-cost raw materials (especially enzymes).
- No standardized platform approach
- DNase/RNase contamination



## Speed to Clinic | Market

- Robust scalable processes
- Standardization of the process
- Process simplification
- Skilled professionals



## Process Improvement


- Increase manufacturing yield
- Reduce cost per dose
- CQA & Analytical tools
- In house capacity expansion



## Flexible Processes

- Multiple product options
- Modular and mobile
- Single use technology
- CDMO service

# mRNA Workflow and Sartorius Solution



## Poll 3

Talk to our experts  
about the mRNA  
related processes  
at

**Raman.Rani@  
Sartorius.com**

# Agenda

- Introduction
- *Upstream – in vitro transcription*
- Downstream Solutions
- Formulation Development & Storage
- Data Analytics



# Unique Challenges of mRNA IVT




Pain points	No standard efficient IVT reaction protocol available	High cost and weak supply chain of raw materials	Lack of dedicated tools supporting PD of mRNA	DNase and RNase contamination / instability of mRNA	No dedicated analytical methods for this new modality
	High-throughput	mRNA expertise	Single-use technology	Analytical solutions	

\*LNP: Lipid nano particle

# In Vitro Transcription Reaction – Production of mRNA from DNA Template

- Reaction time is typically 2-3 hours for batch processes.
- High yield (low reaction volumes)
- IVT is a multi-component reaction :
  - Plasmid (dsDNA)
  - RNA polymerase (e.g. T7)
  - NTPs (optional modified NTPs)
  - Capping reagent (optional)
  - MgCl<sub>2</sub>
  - Pyrophosphatase (optional)
  - RNase inhibitor
  - Spermidine
  - DTT
  - Buffer




Rosa, S. S., et al., Vaccine. 2021 Apr 15; 39(16): 2190–2200.

# Monitoring of IVT: Two Paradigms for Rapid At-Line Analytics


## CIMac Oligo dT - affinity

- One-parameter-at-a-time, faster than Ribogreen
- 'Protein A - mAb' paradigm for mRNA
- Titre of polyadenylated mRNA throughout process





## CIMac PrimaS - multimodal

- Multi-parameter method
- New paradigm for mRNA
- NTP, capping, RNA content. Applies to all RNA modalities









- The IVT reaction can be monitored at-line by **CIMac PrimaS**
- mRNA production kinetics is monitored. Productivity maximum can be identified, to prevent degradation.
- Consumption of nucleotides and concentration of capping reagent can simultaneously be studied
- Effects of feed addition can be studied





Pregeljc, D. et al. "Increasing yield of in vitro transcription reaction with at-line high pressure liquid chromatography monitoring." *Biotechnology and Bioengineering* (2023, 3, 737-747)

# CIMac PrimaS Provides Tight Control Over IVT Reaction



Pregeljc, D. et al. "Increasing yield of in vitro transcription reaction with at-line high pressure liquid chromatography monitoring." *Biotechnology and Bioengineering* (2023, 3, 737-747)



# Optimization of IVT reaction



- Methodological platform for studying individual or combinatorial effects of IVT variables on mRNA production and NTP consumption
- Combined with MODDE® software for DOE and data analysis to derive deeper process understanding

Pregeljc, D. et al. "Increasing yield of in vitro transcription reaction with at-line high pressure liquid chromatography monitoring." *Biotechnology and Bioengineering* (2023, 3, 737-747)

# Batch to Fed-Batch – Monitoring NTPs and mRNA




- PATFix Monitoring of depletion of NTPs; react with feed addition
- Control scale-up of IVT reaction
- Control tech transfers
- Calculate kinetics of NTP consumption → transform to continuous feeding (e.g. AMBR250)

Pregeljc, D. et al. "Increasing yield of in vitro transcription reaction with at-line high pressure liquid chromatography monitoring." *Biotechnology and Bioengineering* (2023, 3, 737-747)

# Combining Fed-Batch With Co-transcriptional Capping to Reduce CoGs


- Combine **fed-batch** and **co-transcriptional capping** paradigms to decrease CoGs by lowering T7 / plasmid requirements
- Feed GTP to keep constant 2 mM concentration
- High capping efficiency (80%) at high mRNA yield (10 g/L) demonstrated with ARCA
- Principle applies to all co-transcriptional cap analogues



Pregeljc, D. et al. "Increasing yield of in vitro transcription reaction with at-line high pressure liquid chromatography monitoring." *Biotechnology and Bioengineering* (2023, 3, 737-747)



# Towards Continuous Production of mRNA




- Determination of optimal IVT conditions and characterisation of reaction kinetics in thermal shaker.
- NTP feeding strategy designed and tested in thermal shaker and transferred to automated scale-up Ambr 250 system.

Skok, J, Megušar, P. et al. "Gram-Scale mRNA Production Using a 250-mL Single-Use Bioreactor." *Chemie Ingenieur Technik* (2022) 94 1928-1935



# Multi-gram mRNA Production in 250 mL Single-Use Bioreactor



- At-line monitoring of NTP in reaction container.
- 2 g mRNA produced from 100 mL starting volume
- cost reduction of up to 50 % per mg of mRNA produced compared to batch mode production by better utilization of enzymes and plasmid.

Skok, J, Megušar, P. et al. "Gram-Scale mRNA Production Using a 250-mL Single-Use Bioreactor." *Chemie Ingenieur Technik* (2022) 94 1928-1935


# Agenda

- Introduction
- Upstream – *in vitro* transcription
- **Downstream Solutions**
- Formulation Development & Storage
- Data Analytics





# Unique Challenges of mRNA Downstream Processing



Pain points	IVT protocol is still evolving, impacting the DSP	There is no established DSP platform that can be applied to all mRNA	Reagents and raw materials are poorly defined and bring contaminants	Impurities are very similar to the target and may interact with it	Many technologies used in labs are not scalable	mRNA are large molecules, shear sensitive and relatively unstable
Needs	Monolithic technology	mRNA expertise	Single-use technology	Proven Scalability	Toolbox of filters	

\*LNP: Lipid nano particle

# What are we Separating mRNA from?

- IVT reaction components (DNA template, enzymes, NTPs, capping and/or other reagents)
- Contaminants in the raw materials:
  - In the plasmid: E.coli proteins, DNA, and RNA.
  - In the enzymes: fragments, host DNA, RNA, proteins.
- Endotoxins, if introduced during processing
- Reagents used for any additional processing, such as plasmid digestion post-synthesis or capping enzymes
- RNA variants: dsRNA, truncations, fragments, aggregates.



mRNA




linear pDNA  
+  
enzymes  
NTPs  
other reagents



by-products

# Chromatography



# Toolbox for mRNA Purification and Analytics



## Advantages of a Toolbox

- Capture mRNA from IVT without need for UF | DF
- Poly-A specific capture with affinity Oligo dT18
- Multimodal capture for RNA without poly-A tail
- Solutions are ready to deploy and screen with new projects or customer wishes.
- Sanitise and re-use CIM® monoliths

# Affinity Capture of Polyadenylated mRNA by CIMmultus Oligo dT18




- IVT > add NaCl > Load on CIMmultus Oligo dT (no concentration, no buffer exchange, no TFF)
- Binding in moderate NaCl concentrations (250 mM – 1.5 M) leads to DBC of 3-4 mg/mL, in Gdn up to 6 mg/mL
- Elution in low concentration buffer or in water

Mencin, N., Štepec, D. et al., Development and scale-up of oligo-dT monolithic chromatographic column for mRNA capture, Separation and Purification Technology, 304 (2023)

Korenč, M. et al, Chromatographic purification with CIMmultus™ Oligo dT increases mRNA stability, Cell & Gene Therapy Insights 2021; 7(9), 1207-1216

# Capture of mRNA Without Poly-A Tail Using PrimaS



- Elution at pH gradient (or step) can separate mRNA from IVT components.
- Capture of mRNA without poly A tail, ability to remove protein, nucleotides, plasmid
- Robust, IVT is applied to the column after initial dilution with loading buffer.




Pyrophosphate  
gradient elution




pH gradient  
elution



# FLD + UV to monitor protein clearance by PrimaS



Pyrophosphatase + mRNA eGFP




T7 polymerase + mRNA eGFP

RNase inhibitor + mRNA eGFP

All proteins + mRNA eGFP

FLD detection demonstrates protein removal in high-salt wash (UV signal too low to detect)  
mRNA elutes in pH gradient

# pH gradient elution from PrimaS does not impact mRNA stability





- mRNA stable for at least 28 days at room temperature after PrimaS capture from IVT

Eluate neutralisation required for mRNA stability.



# Removal of dsRNA by Reverse-Phase Chromatography (SDVB)




- CIMmultus Oligo dT-purified Cas9 mRNA (4000 nt)
- Loaded onto CIMmultus SDVB 8 mL
- Room temperature separation 7.5 → 18% ACN in 50 mM TEAA, pH=7.0
- Removal of dsRNA demonstrated by J2-dot-blot

# mRNA Size Separation Using CIMac SDVB Column

## CIMac SDVB

- Separates ssRNA by size, example below shows RiboRuler High range RNA ladder
- Enables monitoring of shorter transcripts removal, tracking of degradation products during stability studies




CIMac SDVB with 2 $\mu$ m channels, Buffers: MPA (50 mM TEAA + 7.5 % ACN), MPB (50 mM TEAA + 18 % ACN), Method: 0–0.5 min (100 % MPA), 0.5–11 min (gradient to 100 % MPB). Full method not shown. Flow rate 1 mL/min, PATfix™ system, UV absorbance at 260 nm, column thermostat at 60 °C.

# Concentration | Diafiltration | Sterile Filtration

# UFDF for development and manufacturing of mRNA



Sartoflow® Smart TFF System



Ambr® Crossflow  
(5 mL – 500 mL screening)



Select the filters and crossflow devices to get the best recovery from a comprehensive toolbox of scalable solutions:

- Purification to get rid of DNA fragments, residual enzymes such as Polymerase (99 KDa), Dnase (74 kDa), Pyrophosphatase (19 KDa) and other impurities.
- Hydrosart, a state-of-the-art crossflow membrane designed for low adsorption of mRNA and LNP, available in a 50cm<sup>2</sup> and scalable format
- Hollow fiber modules for gentle TFF and maximized mRNA yield in a 52cm<sup>2</sup> scalable format
- Run your experiments using the Sartoflow® Smart benchtop system designed for process development, including a low shear membrane pump for highest mRNA recovery
- Mitigate risk of contamination and RNase degradation with Sartocoon® Self Contained Units, hollow fibers and Maxicaps filters which can be delivered gamma sterile, facilitating closed processing
- Facilitate closed processing by using pre-assembled transfer sets including filters, tubes and connectors



Hydrosart® 100/300  
KD Sartocoon Slice



Hollow Fiber TFF  
modules (up to  
750kD) Green line



Multi-Use TFF  
Systems




Sartocoon Self-  
contained modules  
and assemblies



Single-Use TFF Systems

# UFDF for mRNA and LNP



Parameters	Value
Device	Hollow fiber   Hydrosart (100   300 kDa)
TMP	< 0.5 bar
Trial loading (L/sqm)	Up to 70 L/m <sup>2</sup>
Volumetric concentration factor	8X
Diafiltration volume	Up to 12 DV

Hydrosart 300 kDa performance data for LNP

- For LNP no change in Particle size were observed
- Cassettes are tested up to 10 kb size mRNA

# Sterile Filtration for mRNA and LNPs



Screening

Sartoscale 25 device

PES

Sartopore<sup>®</sup> 2/XLG/Platinum



A toolbox of filters to find the best balance between filtration capacity and recovery

- Broadest range of PES membrane for high mRNA and LNP recovery
- PES filters performance better than other MOC
- Low non-specific binding which results in high recovery
- Simplify process development with the use of scalable options
- Off-the-shelf standardized pre-designed solutions
- Ensure process automation and control by using single-use systems



SU Assemblies

Filter transfer sets & Maxicap MR

## Poll 4

Talk to our experts  
about the mRNA  
related processes  
at

**Raman.Rani@  
Sartorius.com**

## Poll 5

Talk to our experts  
about the mRNA  
related processes  
at

**Raman.Rani@  
Sartorius.com**



## Poll 6

Talk to our experts  
about the mRNA  
related processes  
at

**Raman.Rani@  
Sartorius.com**

# Agenda

- Introduction
- Upstream – *in vitro* transcription
- Downstream Solutions
- **Formulation Development & Storage**
- Data Analytics




# Unique Challenges of LNP Development, Storage and Shipping of mRNA



Pain points	LNPs are fragile and shear sensitive	mRNA is relatively unstable and easily degraded by RNase	Performance differs from small-scale to late-stage	Current LNPs are stored at extremely low temperatures
Needs	Toolbox of filters	Single-use technology	Proven Scalability	Controlled Freeze Thaw

\*LNP: Lipid nano particle


# Delivering via Lipid Nano Particles (LNP)



- Lipids composition of LNP may vary to optimize cellular uptake, endosomal escape, mRNA payload
- LNP size varies from 50 – 150 nm

## Challenges

- Aggregation and Stability
- Control of lipid type, source and quality



## The typical liposome for a mRNA vaccine contains 4 lipids

- 1x Cationic lipid
- 1x Helper lipid
- 1x Cholesterol
- 1x PEG lipid

## Cationic/ionizable lipids

e.g. DOTMA/DOTAP or proprietary lipid

- Nucleic acid complexation
- Membrane fusion

## “Stealth” PEG lipids

- PEG 2000
- Hydrophilic surface
- Steric hindrance

## Structural helper lipids

e.g. DSPC, DPPC

- Bilayer support

## Cholesterol

- Integrity
- Endosomal release

Reference : Rein Verbeke, et.al , Three decades of messenger RNA vaccine development, August 2019, [Nano Today](#) 28(Pt 1):100766

# Simplify Formulation and Filling Development of mRNA

## Speed up your LNP development with high throughput screening TFF and controlled Freeze/Thaw system

- ✓ Screening and identification of CPPs and CQAs
- ✓ Set Design Space

### LNP formulation and filling development mRNA package

#### TFF systems & consumables:

- [Sartoflow® Smart TFF System](#)
- [Hydrosart® 100/300KD](#)  
[Sartocon® slice 50 TFF Filters](#)
- [Hollow Fiber TFF](#)  
[\(100KD/300KD\)](#)

#### Filters:


- [Sartopore® 2/XLG Sterile Filters](#)

#### Software:

- [MODDE® DOE Software](#)

#### Freeze/Thaw System:

- [Celsius® S<sup>3</sup> Benchtop Freezer](#)



# Controlled Freeze/Thaw Characterization for mRNA Formulation



Celsius® S<sup>3</sup> Benchtop Freezer

Celsius® Pak Bags



Ensure stability of drug product and drug substance during storage by optimizing the freeze & thaw steps:



- Benefit from an end-to-end integrated approach for faster and more efficient testing and validation of freeze & thaw, assuring quality at scale
- Investigate optimal process conditions using the Celsius S3, generate a consistent samples library to investigate product stability, storage and shipping process steps
- Enables freeze/thaw characterization in 30–100 mL single-use bags with the same product contact material from lab-scale to commercial-scale Celsius® platforms
- Automated operation and data collection, combined with accurate temperature monitoring, at well-defined last-point-to-freeze location

# Late-Stage mRNA Storage and Shipment with Controlled Freeze|Thaw

## Celsius® CFT Freeze and Thaw Solution



Celsius® FT Plate Freezer



Celsius® Pak Bags



Celsius® FS Filling Station



Celsius® SSM Shipper



Celsius® FFTp SU Container



Biosealer® TC Sterile Disconnection



Biowelder® TC Sterile Connection



Ensure safe and simple preservation and transfer of frozen mRNA and LNPs:



- Control your process to minimize the adverse effects of the cryoconcentration in your LNP with a controlled-rate freeze and thaw system with Celsius® CFT
- Ensure consistent freeze|thaw performances in single-use bags with same product contact material from 1L up to 100L. One film for all to facilitate process validation with Celsius® CFT
- Celsius® FFTp - pre-assembled and ready to fill single-use containers to leverage existing infrastructure, fully qualified down to -80°C
- Mitigate the risk of contamination and RNase degradation with gamma-irradiated single-use bags
- Sterility Assurance during connection and disconnection, thanks to effective solutions adapted to your process requirements
- Confidence® Validation Services to support your DNase|RNase mitigation strategy.


# Agenda

- Introduction
- Upstream – *in vitro* transcription
- Downstream Solutions
- Formulation Development & Storage
- **Data Analytics**





# Maximize Process Understanding & Control Using Data Analytics



# Sartorius Solutions Simplify Progress And Fast-Track Time-to-Market

## Expertise




Combining high-throughput with innovative analytics and supportive experts enable you to quickly develop in-house expertise

## Innovation



Only innovative solutions dedicated to mRNA applications can eliminate some of your challenges and support your next generation process

## Scalability



Develop the process with the end in mind to ensure successful implementation at late-stage from a cost, process and supply perspective

# Thank you!

**Niddhivinayak S Mishra**  
FAS Manager, BPS India  
Sartorius Stedim India Pvt Ltd.  
Headoffice – Gottingen, Germany  
Mobil +91 8884001204  
Niddhivinayak.Mishra@sartorius.com

**Sekirnik, Rok**  
Head of Profit Center  
BIA Management  
Sartorius Stedim  
Headoffice – Gottingen, Germany  
rok.sekirnik@sartorius.com

The Sartorius logo is displayed in a bold, black, sans-serif font against a bright yellow background. The letters are closely spaced and have a clean, modern appearance.