

Centrifuges in vaccine production processes

GEA Separation technology

THOMAS MÜLLER, 06.08.2020, OELDE, GERMANY





Thomas Mueller

Global Product Manager Sales
Pharma Biotechnology

With GEA since 2006
specialized in separation technology in pharma processes

Formerly 15 years with pharma company Baxter, starting as automation specialist, later taking on various positions as pharma specialist



- 50 years old, married, two sons aged 16 and 19
- Main education: automation technician
- My second life is to work around lifesaving (Baywatch)

Contact: Thomas Müller, Tel.: +49 1721615934, thomas.mueller@gea.com

- About GEA
- Separation Technology
- Biosafety
- Separation technology for vaccine processes
- Experience with vaccine productions

In touch with GEA every day - Our applications

Dairy Farming and Processing



Approx. one
quarter of
processed milk
comes from
GEA production
systems

Food



Every third
chicken
nugget is
produced
using GEA
technology

Beverages



Approx. every
second liter
of beer is
brewed with the
aid of systems
and process
solutions from
GEA

Pharma



Every fourth
liter of human
blood for
making
plasma-derived
products is
processed
using GEA
equipment

Chemical



More than one
third of all
polymer
producers
are using
GEA drying
technology

Utilities



Each industry
we serve
utilizes
industrial
refrigeration
technology
from GEA

Marine



Every second
container ship
in the world
sails with
GEA marine
equipment on
board

GEA is one of the largest suppliers of process technology to the food industry and to a wide range of other industries.



4,751

million order intake (EUR)



17,863

employees (FTEs)

The international technology group focuses on process technology, components and sustainable energy solutions for sophisticated production processes in diverse end-user markets.




564

million operating
EBITDA (EUR)



4,605

million revenue (EUR)



1.31

earnings per share (EUR)


The company is listed on the German MDAX stock index (G1A, WKN 660 200) and included in the STOXX® Europe 600 Index. In addition, the company is listed in selected MSCI Global Sustainability Indexes.



12.2

% operating EBITDA margin

Strong local presence around the globe




GEA centrifuges

GEA
SEPARATION
126
YEARS
Since 1893

The first **centrifuge** was built in Oelde in **1893** - and we've been at the forefront of mechanical separation technology ever since.

The **separators**, **decanter**s and **process lines** prove their value worldwide in high-performance companies of all shapes and sizes, sometimes under extreme operational conditions.



Centrifugal separation: approx. 3,500 applications



Chemical/Pharma



Renewable Ressources



Beverage



Dairy



Environmental



Oil & Gas, Energy



Marine



Centrifugal Separators

- With bowl diameters of 160 mm to 1,050 mm
- Throughput up to 500,000 l/h
- Particle size from 0.5 μm
- “Liquids oriented”



Decanter Centrifuges

- With bowl diameters of 200 mm to 1,030 mm
- Throughput up to 350,000 l/h
- Particle size from 5 μm
- “Solids oriented”



**Pharmaceutical
skids perfectly
designed for
standard and
customized
requirements**

Maximum safety,
maximum yield.

Fermentation



Technologies

- Fermenters
- Bioreactors

Liquid dosage applications

- Vaccines
- Insulin
- Growth Hormones
- mAbs
- Therapeutical Proteins

Homogenization




Technologies

- Cell Disruption
- Homogenization

Liquid dosage applications

- Vaccines
- Therapeutics
- Diagnostics
- Enzymatic Preparations

Separation



Technologies

- Mechanical Separation

Liquid dosage applications

- Vaccines
- Hormones
- Human blood proteins
- Mammalian cell culture
- Starter cultures
- Enzymes

Safety



Efficiency




Biosafety is of highest
priority



Safety: Strictly aseptic process

Absolutely reliable steam sterilization required for vaccine and sera:

- Compliance with FDA requirements
- Compliance with GMP requirements
- Sterile complete closed system
- Guarantee on the employed sterile separator skids
- CIP/SIP including biocontainment





Securely clean with CIP


All GEA Pharma Separators

Absolutely clean without dismantling


- **Fully automatic CIP system available**
The cleaning cycle is controlled by a programmable logic controller. The program sequence for cleaning can be adapted according to local requirements.

Securely sterile design with steam SIP


- **A sterile, complete closed system ensured**
according to the FDA specifications
(Food and Drug Administration)
- **Proven**
in more than 800 SIP Skid Units installed world-wide




Our separation solutions for vaccine processes



Comparison: “Live Vaccine” und “Non-live Vaccine”




Live Vaccine



Non-live Vaccine

Vaccine Process – Human vaccines



Human vaccines

Cell separation

- 1 Microorganism / Substrate / Aeration
- 2 Fermentation
- 3 Biomass

Non-Live Vaccine Processing

- 4 Harvest tank
- 5 Kill tank
- 6 Further processing

Live Vaccine Processing


- 7 Live vaccine to freeze dryer
- 8 Additives

a Self-ejecting or nozzle separator


- 9 Waste

- Yellow arrow: supernatant
- Brown arrow: solids

Vaccine Process – Veterinary Vaccines (FMD)



Selfdesludging centrifuge



Requirements – GEA features – Your benefits

Main features	Your main benefits
Patented hydrohermetic product feed inlet	Gentle product treatment tested on shear sensitive mammalian cell cultures for best product quality
Unique patented gap-free welded disc spacers	Highest hygienic standard
Variable ejection volume	Adaption to different product properties
Discharge of separated solids with the patented GEA hydrostop system	Precise, reproducible and step less adjustable ejection volumes

Main features	Your main benefits
Newest generation for HMI and PLC series by SIEMENS	Long availability of spare parts
Combined compact control cabinet for MCC and PLC mounted on the skid	Pre-cabled and pre-tested
GAMP 5*	Guideline for programming of control system for pharmaceutical industry including the relevant documentation

Separator



Automation



Requirements – GEA features – Your benefits


Main features	Your main benefits
Complete skid including piping and control cabinet	Plug-and-play
Customized footprint and package assembly	Adaptation to customers' needs
Pretesting of the complete skid on our pharmaceutical test bay with purified water and clean steam	Reduces time for testing during commissioning

Main features	Your main benefits
Qualification documents of complete scope of supply done by GEA	Prequalified system, easy to implement in customers' qualification files
GEA verification plan including IQ/OQ documents filled out and tested	Minimized qualification work on customers' side

Package unit



Qualification



Special requirements for documentation

- Legal requirements for the pharmaceutical industry



- Comprehensive requirements of the pharmaceutical customers

- Machines, equipment, spare bowls and spare parts



FAT - SAT

Product recipes

ASME BPE

21 CFR part 11

X-Ray



Filter integrity test

GAMP 5 - HDS - SDS - FDS

DQ - IQ - OQ

HAZOP Meeting

Monthly meetings / telcons

Riboflavin test acc. GEA standard

Centrifuge Product Range

Self cleaning separators in steam sterilisable design GEA Pharma separator aseptic

Machine type	Product capacity depending on product and process, in l/h
CSC 6	100 - 200
CSC 15	250 - 500
CSC 20	500 - 1.000
aseptic 30*	500 – 2.500
aseptic 100*	1.800 - 3.700
aseptic 170*	3.000 - 6.000
aseptic 300*	6.000 - 12.000

* ASME BPE design and Integrated direct drive (IDD)



Efficiency: Skid installation for short ROI

- Skid units ready for plug and play
- Pilot skids to set up new processes (option of easy up scaling)
- Large variety of applications
- Large variety of product variants
- Sophisticated service along the life cycle
- Interdisciplinary team of experts supporting GEA project manager



Experience with vaccine productions



Sold 160 package units



Until
2019

50 customers

CSC 6 Package Unit with Automatic SIP



CSC 15 Package Unit with automatic SIP for the Recovery of Vaccines




Sold 116 package units




Until
2019

38 customers

CSC 6 Package Unit for the Recovery of Veterinary Vaccines



CSE 80 Package Unit for the Recovery of Vaccines (typical for FMD)



- Global technology leader
- Specialist for pharma applications
- Experts on safety issues
- Experts on efficiency





engineering for
a better world