

Energizing pharma

Efficient and secure: Vaccine Production with GEA pharma separators

DR. MARKUS KÜHBERGER

SENIOR DIRECTOR BUSINESS LINE CHEMICAL & PHARMA

GEA at a **Glance**

"GEA is one of the world's largest systems suppliers for the food, beverage and pharmaceutical sectors. A major focus is on continuously enhancing the sustainability and efficiency of customers' production processes. GEA plants, processes and components help achieve significant reductions in carbon emissions, plastic use and food waste in production worldwide. In this way, GEA makes a decisive contribution toward a sustainable future, fully in line with its corporate philosophy of "engineering for a better world."

GEA is listed on the German MDAX and the STOXX[®] Europe 600 Index and also included in the DAX 50 ESG and MSCI Global Sustainability indexes.



Local presence worldwide





Separation solutions in human vaccines



Non-Live vaccines

- Influenza
- Cholera
- Bubonic plague
- Hepatitis A or hepatitis B
- Whooping cough (pertussis)
- Tetanus
- Diphtheria
- Pneumococcal bug

Culture media \Rightarrow Valuable phase

GE?

engineering for a better world



Live vaccines

- Measles
- MumpsRubella
- Nubella
- Yellow fever
- Chickenpox
- BCG (= tuberculosis vaccine)
- Typhus
- Human rotaviruses
- Whooping cough (pertussis)

Culture media ⇒ Waste



The new GEA Pharma separator lines aseptic & pure



GEA Pharma Separator aseptic



GEA Pharma Separator pure



2 full separator lines	Multifunctional machine by GEA flexChange concept
ASME BPE compliant design	Optimized Bowl concept including the
Utmost variability (flexibility)	NEW GEA flexicon nozzle separator
Integrated Direct Drive concept	Skid configurator

Selection process (Customer requirements)



Cleanability, surface texture and documentation always adaptable to your needs



Let's get into technical details





Maximum yield thanks to high precision desluding



GEA self-desludging centrifuges



Typical feed conditions in vaccine productions:

3-5 % solids in volume, depending on vaccine30 -1000 l/h, up to 7000 l/h possible for larger machines





Gentle product handling, higher viability of cells, more efficient clearification

Hydrohermetic Feed patented by GEA Westfalia Separator Group

Common Feed Product accelerated freely Hydrohermetic Feed -

Product accelerated within the liquid, underneath the disk



No turbulences, no splashes, no shear forces

Cleanability – CIP and SIP



Total bowl ejections -

possible with hydrohermetic feed system





Spray Nozzles



Steam Sterilization

The next level of pharmaceutical cleanability: Patented gap-free welded disc spacers

GEA's **patented laser-welded** disk spacers lead you to the next level of pharmaceutical cleanability.

Perfectly clean: High quality laser-welded spacers Maximum welding precision leaves no dead spaces in the disk stack, providing the basis for a highly effective CIP and SIP cleaning.

Optionally available for GEA pharma separator pure.





Unique: GEA Pharma test bay for pre-qualification





Pre-qualification of separator centrifuge skids

Flexibly and safely on the GEA pharma test bay.

Reference Customers Human Vaccines



Sold 172 package units









Single-use GEA kytero[®] Separator

All our competence – in an easy-to-operate single use skid





GEA kytero[®] Separator: 1 new design – 3 × benefits for you

Easy and time-saving operating concept

Maximum-purity GEA breezeDrive[®] revolution

Proven high-yield disk stack technology





I am happy to answer your questions:





Dr. Markus Kühberger

Senior Director Business Line Chemical & Pharma GEA Westfalia Separator Group GmbH Werner-Habig-Str. 1 D-59302 Oelde /GERMANY

 Email
 markus.kuehberger@gea.com

 Mobile
 +49 157 805 98 211

 Landline
 +492522775545

 Web
 www.gea.com