

Webinar: Polio Analytical Services



Agenda

- Introduction Viroclinics-DDL
- Background Polio
- GAP III Regulations
- Analytical assays in a Polio High Containment Facility
- Q & A



Introduction Viroclinics-DDL



European HQ and Laboratory Sites



Viroclinics

Clinical diagnostics laboratory.
Support staff and corporate HQ
location.



Viroclinics Xplore

Early stage and preclinical
services. Capabilities at BSL2
and BSL3(+) level.



DDL Diagnostic Laboratory


Focused on molecular diagnostics,
molecular biology and molecular
pathology.



Bridged Assays

Based in Shanghai, central virology
lab service for phase II, III and IV
clinical trials.

Our Therapeutic Fields



Infectious Diseases

Respiratory viruses


Influenza
RSV
MERS-CoV
SARS-CoV
SARS-CoV-2 (COVID-19)

Blood-borne viruses

Hepatitis B, C, D
HIV

Other viral targets

HPV, **Polio**, Dengue,
Chikungunya ... and more




Gene Therapy

Adeno-associated viruses

AAV-2
AAV-9
... and more



Bacteriology



Immuno-Oncology

Under development



Microbiome

BSL 3 Facilities

CUTTING
EDGE LAB

BSL 2

BSL 3

State-of-the-Art
Facilities




DEDICATED
SCIENTISTS &
CONSULTANTS



In-house preclinical and clinical laboratories

Our Capabilities - Preclinical Services

- 1600 m²
- SPF, Conventional, BSL2 and BSL3(+)
- GAP III Polio Certificate of Participation
- Negatively-pressurized isolator cages
- Range of different species
- Solutions for any kind of experimental set up: challenge and transmission studies





Our Capabilities - Preclinical Services

Type of research


- Vaccine Efficacy
- Passive Immunization
- Immunogenicity
- Pharmacokinetics
- Immunopathology
- Pathogenicity
- Antiviral Therapy
- Antiviral Prophylaxis
- Safety
- Drug Delivery
- Resistance / Mutants
- Transmission
- ...*(more)*

Read Outs

- Virus Titers
- Virus Copies
- Antibody
- CMI
- Genomics
- Pathology
- Histopathology
- Mortality
- Clinical Signs
- Imaging
- ...*(more)*



Background Polio




Poliovirus

Picornaviridae

(+) ss RNA genome

~7500 nucleotides

30 nm



3 serotypes:

Poliovirus 1 (PV1)

Poliovirus 2 (PV2)

Poliovirus 3 (PV3)

Wild poliovirus (WPV)

Circulating Vaccine derived
poliovirus (cVDPV)

Image: www.virology.ws

polio | **global**
eradication
initiative

Polio vaccines

Inactivated polio vaccine (IPV)

- Jonas Salk, 1955
- Contains inactivated viruses of all 3 serotypes
- Prevents disease
- More expensive
- Injectable, requires skilled personnel



The "Biltoven Unit," designed by the Dutch microbiologists van Hemert (pictured) and van Wezel, generated large quantities of poliovirus for vaccine production in the 1960s.

Oral polio vaccine (OPV)

- Albert Sabin, 1963
- Contains attenuated viruses of all 3 serotypes
- Prevents disease
- Cheaper
- Easily administered, passive immunization



Oral polio vaccination with attenuated virus



**Vaccination of
naïve individual**



replication in the intestine




**Immunized
individual**

Oral polio vaccination with attenuated virus



**Vaccinated
individual**



**Excretion and spread to others
in close contact**




**Passive
immunization**


Circulating vaccine-derived poliovirus



community with
insufficient vaccine
coverage




Prolonged circulation and
spread, reacquired
neurovirulence



Vaccine-derived
poliovirus
outbreak

Poliomyelitis Global annual reported cases



Polio eradication

2019

Two out of three wild poliovirus strains eradicated

Global eradication of wild poliovirus type 3 declared on World Polio Day

24 October 2019


In an historic announcement on World Polio Day, an independent commission of experts concluded that wild poliovirus type 3 (WPV3) has been eradicated worldwide. Following the eradication of smallpox and wild poliovirus type 2, this news represents a historic achievement for humanity.

2020

Africa Kicks Out Wild Polio

On August 25, 2020, the African Region was officially certified as wild poliovirus-free. The Centers for Disease Control and Prevention (CDC) has served as the lead implementing partner for U.S. efforts to eradicate polio, playing a pivotal role for three decades in helping African countries and the continent reach this milestone. The goal of our work has been to protect all children from the dangerous, sometimes deadly, yet entirely vaccine-preventable consequences of polio transmission and outbreaks.

Risks and challenges for polio eradication



**Spread of wild poliovirus
from endemic countries**

- WPV1 cases
- cVDPV2 cases

Risks and challenges for polio eradication



Vaccine derived poliovirus outbreaks

- Lack of basic health infrastructure that limits vaccine distribution
- Neglect of polio vaccination programs

- cVDPV1 cases
- cVDPV2 cases

Risks and challenges for polio eradication

Laboratory escapes

Polio virus escaped twice from laboratories

TNN | Oct 6, 2004, 01:25 IST

Belgium: Live polio virus solution accidentally released into local water

by ROBERT HERRIMAN

🕒 September 12, 2014

📍 Europe

💬 19 Comments

HEALTH

POLIO VIRUS FOUND IN SEWAGE AT PUBLIC HEALTH INSTITUTE

By Janene Pieters on Wednesday, 2 September 2020 - 11:50




Employee vaccine manufacturer infected with polio virus

Publication date 05/02/2017 - 00:00

GAP III

WHO Global Action Plan to minimize poliovirus facility-associated risk after type-specific eradication of wild polioviruses and sequential cessation of oral polio vaccine use.

Poliovirus type 2, type 3 materials



Data reported by WHO Regional Offices as of 22 January 2019 and subject to change


The boundaries and names shown and the designations used on this map do not imply the expression of any opinion whatsoever on the part of the World Health Organization concerning the legal status of any country, territory, city or area or of its authorities, or concerning the delimitation of its frontiers or boundaries. Dotted and dashed lines on maps represent approximate border lines for which there may not yet be full agreement.

 World Health Organization
© WHO 2016. All rights reserved.

26 countries plan to retain poliovirus type 2/3 materials in 78 designated poliovirus-essential facilities (PEFs)


Viroclinics Xplore has established a GAP III Polio Certificate of Participation

Containment




Adapted: World Health Organization 2019

BioRisk Management



Adapted: World Health Organization 2019

Analytical assays in a Polio High Containment Facility




Assays


- Virus neutralization
- Quantitative virus detection
- Preclinical testing
- Neurovirulence testing
- Deep sequencing

Virus neutralization

Assay performed in GAPIII laboratory:




Advantage over Crystal Violet readout:



Titration in cell culture and preparation of inocula


1. Dilution of samples

1 log₁₀ and 0.5 log₁₀ dilution steps




2. Transfer of dilutions

To 96-well plate




Crystal violet staining



3. Addition of Hep2 cells

To 96-well plate



WST-8



4. Incubate

- At 35°C, without CO₂
- Using sealing films for 96-well plates
- Incubate for 6 days

Outlined in the WHO Recommendations to assure the Quality, Safety and Efficacy of Live Attenuated Poliomyelitis Vaccine

Mouse Model for Neurovirulence Testing



TgPVR21 mouse line



CENTRAL INSTITUTE FOR EXPERIMENTAL ANIMALS
3-25-12 Tonomachi, Kawasaki-ku, Kawasaki, 210-0821 JAPAN
Phone: +81-44-201-8510 Fax: +81-44-201-8511

Certificate

Viroclinics Biosciences
Marconistraat 16
3029 AK Rotterdam
The Netherlands

Date : 26th September, 2016

Central Institute of Experimental Animals (CIEA) for WHO has an obligation to provide suggestions for appropriate maintenance and use of TgPVR mice to prevent the escape of TgPVR mice as a possible Porio virus carrier. Especially, for a new facility (country) to start maintenance of TgPVR mice, CIEA must check a sketch map and pictures of animal facility and Standard Operating Procedure (SOP) before using TgPVR mice.

According to the World Health Organization (WHO) guideline reported in November, 1992, the inspection of your animal facility (Lander Campus: Nistelrooise baan 3, 5374 RE Schaijk, The Netherlands) has been performed. For the inspection, the following 4 major points of WHO guideline, 1) Animal facility, 2) Rearing cage, 3) Rack, 4) Care and management of TgPVR mice were examined.

As the result of the inspection, we approved that your animal facility was appropriate for maintenance and use of TgPVR mice.


Sincerely yours,



Ryuta Nomura,
Executive Director
Central Institute for Experimental Animals

Neurovirulence Test

Needle inserted between the spinous process and the first lumbar vertebrae




Stages	Physical signs		
Normal	Grips the edge of the cage	Walks normally on the grid and on a flat surface	Full ability to move limbs forward
Weak	Unable to grip the edge of the cage	Walks normally on the grid and on a flat surface	Full ability to move limbs forward
Paresis/ Partial paralysis	Unable to grip the edge of the cage	Limbs fall through the grid more than once while walking and toes curl repeatedly while walking on a flat surface	At least a partial to move limb forward
Paralysis	Unable to grip the edge of the cage	No use of limb on grid or flat surface	Inability to move the limb forward

Intraspinal Injection Training with Indian Ink



Poliovirus Genome

- (+) ssRNA virus
- ~7500 nt genome




Quasi-Species

Polymerase (RdRP, 3D^{pol}) lacks proofreading activity


- Mutations
- Recombination

→ Evolution
→ Pathogenesis
→ Escape mutations /
treatment resistance



Determine sequence of virus (quasi-species) in a sample


Sequencing Techniques: 1. Sanger Sequencing



- Determine consensus sequence
→ average of all genomes in a pool
- (sensitive and quantitative)
determination of variants or mutations
not possible

Sequencing Techniques: 2. Next-Generation Sequencing (Illumina Technology)

- Sequence **thousands** of genomes **individually** per run
→ Deep sequencing
- Most common base at each position
→ Consensus sequence
- Differences between individual genomes
→ Variant




Poliovirus Next-Generation Sequencing at Viroclinics

- Samples
 - Clinical samples (e.g. stool)
 - Cell culture isolates
 - Vaccine lots
- Lysis → in GAPIII facility
- RNA isolation, cDNA synthesis & whole genome amplification
- Library preparation
- Illumina sequencing
- Data analysis using in-house pipeline




Results

- Mapping to reference (single/multiple references)
- Variant calling




Results (Complex Samples)

Mixed samples:



Results (Complex Samples)

Recombinants between types:



Our Track Record

OUR TRACK RECORD

700
Studies



OUR TRACK RECORD

Within the last 3 years we processed

150,000
Samples



OUR CONTRIBUTION TO SCIENCE

More than **40**
Publications in peer-
reviewed journals



OUR TRACK RECORD

100%
On Time
Project Delivery



GLOBAL OPERATIONS

INTERNATIONAL
Presence



OUR PERFORMANCE

Within the last 2 years we processed

15,000
Phenotypic Assays



OUR PERFORMANCE

Within the last 2 years we processed

10,000
Genotypic Assays



Post Study Sample Storage



Sample integrity guaranteed with back up systems and 24/7 security and remote monitoring

Polio High Containment Facility

- Virus titration/quantification using fully validated CCID50 protocols
- Neurovirulence testing through a certified and validated assay in transgenic mice (TgPVR21)
- Virus neutralization fully validated on 6 strains
- Transgenic mouse model available for preclinical testing
- Mutational analysis of (whole genome) poliovirus through Next Generation Sequencing
- Top class logistical services to get the samples to our facility
- Storage of samples and (vaccine) virus stocks
- Research and Development



18,000 ft² Preclinical Laboratory Area

Viroclinics contributes to polio eradication by offering its services to Vaccine manufacturers

Questions?



Contact:

Ivo Ploemen
Consultant BD
ploemen@viroclinics.com

For additional info, please visit:
www.viroclinics.com



Viroclinics
Rotterdam Science Tower
Marconistraat 16
3029 AK Rotterdam
The Netherlands
www.viroclinics.com

DDL Diagnostic Laboratory

Visseringlaan 25
2288 ER Rijswijk
The Netherlands
www.ddl.nl