



# Vaccine Processing – an overview

Dr. Mats Lundgren


GE Healthcare Life Sciences

Imagination at work


# How Vaccines are manufactured

## The Vaccines


Bacteria based




Virus based



Protein based



Polysaccharide based



DNA based



## The Manufacturing process

Cell culture / Fermentation



Purification




Fill and Finish


Analysis (QC/QA)

Number and order of the different steps depends on the specific vaccine production

# Different types of marketed influenza vaccines




Whole virus



Split virus



Subunit



Live attenuated

# The evolution of vaccine processes

## 1st generation processes:

*Focus on upstream, optional inactivation*

## 2nd generation processes:

*Separations based on centrifugation, filtration*

## Currently developed processes:

*Quality based approach: Quality by Design*

*Focus on process understanding of entire process incl. purification and virus safety*



# Outline of presentation

Cell substrates for virus production

Cell culture using Microcarriers

Downstream purification of vaccines

Modernizing legacy Vaccine processes

Conclusions




# Cell substrates for virus production



# Selecting a cell line for virus production

- Cell substrate evolution from primary to diploid to continuous cell lines...
- Modern options: Vero, MDCK, EBx™, AGE, PER.C6™ ...
- Requirements
  - Suitable for GMP production
  - Good safety track record
  - Good virus propagation
  - Broadly and highly permissive
  - Scalable to high volume production



from: Pereira et al. Biotech Bioeng; 2004; 85; 5



# MDCK and Vero cells

	MDCK	Vero
+	<ul style="list-style-type: none"><li>• Higher productivity</li><li>• Technically easier</li><li>• Less risk for propagation of adventitious viruses</li></ul>	<ul style="list-style-type: none"><li>• Platform cell line (can be used for several virus vaccines)</li><li>• Good safety record</li><li>• Used for several marketed vaccines</li></ul>
-	<ul style="list-style-type: none"><li>• Potential tumorigenicity/ oncogenicity</li><li>• New cell substrate</li><li>• Restricted to influenza</li></ul>	<ul style="list-style-type: none"><li>• Lower productivity</li><li>• Technically challenging</li><li>• Potential propagation of adventitious viruses</li></ul>





# Virus safety

## EP citation:

*“Seed lots/cell banks. The master seed lot or cell bank is identified by historical records that include information on its origin and subsequent manipulation. Suitable measures are taken to ensure that no extraneous agent or undesirable substance is present in a master or working seed lot or a cell bank.”*



# Cell culture using Microcarriers



# Scale up of adherent cell cultures

Increase volume



Increase number of units



Genetic Engineering News, 2007

One 2500 L bioreactor with a carrier concentration of 3 g/L (Cytodex™ 1) provides the same surface area as 40 000 roller bottles (850 cm<sup>2</sup>/bottle)



# Viruses produced in microcarrier cultures

Adenovirus

Bovine rhinotracheitis

Endogenous C type

Equine rhinopneumonitis

Foot and mouth

Group B arboviruses

HAV

Herpes

Influenza

Japanese encephalitis

Marek's

Papova virus

Polio

Polyoma

Pseudorabies

Rabies

RSV

Rous sarcoma

Rubella

Sendai

SV40

Sindbis

Small pox

Vaccinia

Vesicular stomatitis



# Cell culture media and serum

Serum - Ensure quality, traceability and origin

Classical media


Animal origin free media

Complex media containing hydrolysates


Chemically defined media




# The effect of cell culture media



Medium 1



Medium 2



Medium 3



# Bioreactors – Fixed vs Disposable

## *Control and scalability*



Stainless  
steel



WAVE



10L

2000L

XDR





20 8 2003





# Large scale vaccine production Baxter Biosciences

EC GMP licensed BSL3 (Sept 2004)  
20 million doses plant  
Vero cells on Cytodex™ in protein free  
medium – 6000L scale

Presented at the conference „Influenza Vaccines for the  
world“, Vienna 2006

20 8 2003



# Downstream purification of vaccines




# Downstream processing of viruses

## Available technologies

### Harvest

- Lytic virus
- Non-lytic virus
  - Detergent
  - Mechanical disruption / Homogenization
  - Osmotic shock
  - Freeze-thaw




# Safety and quality is priority

## Regulatory requirements


- Safe vaccine with no or minimal adverse effects
- Effective dose
- Stability
- Process control
- Reproducible process



# Impurity challenges after lysis



# Goal with purification




# Downstream processing of viruses

## Available technologies

### Clarification

- Filtration
  - Normal flow
  - Tangential flow
- Centrifugation




# Downstream processing of viruses

## Available technologies

### Primary purification

- Tangential flow filtration (TFF)
- Density gradient centrifugation
- Precipitation
- Chromatography






# Downstream processing of viruses

## Available technologies

### Secondary purification

- Density gradient centrifugation
- Selective precipitation
- Chromatography
  - IEX, MM, AC, HIC, SEC
  - Bead format (Packed bed)
  - Membrane format (Capsule)





# Ion exchange chromatography

## Anion exchange chromatography


- (-) Negatively charged molecules binds to (+) positively charged ligands

## Cation exchange chromatography

- (+) Positively charged molecules binds to (-) negatively charged ligands



# Size exclusion chromatography




Excluded from pores



Enter a fraction of the pores



Enter all pores



# Affinity chromatography


## Specific binding

Few affinity resins available for vaccines


- Agarose based affinity resin for adeno associated virus
- Pseudo affinity resins for influenza
  - sulphated cellulose
  - sulphated dextrane



# Chromatographic purification of large molecules can be challenging




# Core bead chromatography



- Host cell proteins and DNA fragments bind to the core and viruses stay in the void.


# Process example

## Polio IPV



# The history of Polio vaccine processes

- 1955: Inactivated Polio vaccine (IPV) launched (Salk Type)
- 1960: Attenuated Polio vaccine launched (Sabin type)
- 1960s: Collaboration between Prof. Van Wezel (RIVM/NVI Netherlands) and GE (former Pharmacia) around microcarrier cultures of primary monkey cells.
- 1970s: New IPV purification method using chromatography resins
- 1980s: Switch to Vero cell production
- 2010s: Updating the IPV processes using modern technology





# Modernizing legacy Vaccine processes



# A Modern Solution for Acellular Pertussis Vaccine

The background of the slide is a dense field of green, rod-shaped bacteria, likely representing pertussis bacteria. The bacteria are oriented in various directions, some appearing in sharp focus while others are blurred in the background, creating a sense of depth. The overall color palette is a range of green tones, from light to dark.

# Whole-cell (wP) - Acellular Pertussis (aP)

## wP Vaccines

70 year old technology based on killed *B. pertussis* strains

High protection efficiency ~78%

Associated with side effects and safety concerns

The reactogenicity of wP vaccine was thought to be too high to permit routine use in older children, adolescents and adults.

## aP Vaccines

Introduced in 1990's

aP contain  $\geq 1$  of the separately purified antigens: pertussis toxin (PT), filamentous hemagglutinin (FHA), pertactin (PRN), and fimbriae (FIM) type 2 and 3.

aP is now the dominant type in the industrialized world

aP containing vaccines with reduced concentrations of the antigen have been formulated for use in adolescents and adults



# Project Goal

## Traditional process

Chinese pharmacopeia requirement and current situation

- Contain 2 antigens:  
Pertussis toxoid (PT),  
Filamentous Hemagglutinin (FHA)
- Purity >85% (SDS-PAGE)

Yield around 10%

Lack of stable antigen quantitative assay

## Current Project

Develop a modern process for pertussis vaccine

- Contain 3 antigens:  
Pertussis toxoid (PT),  
Filamentous Hemagglutinin (FHA )  
Pertactin (PRN)
- Purity >95% (SDS-PAGE)

Yield >30%

Establish quantitative antigen determination using Biacore™ platform




# Traditional Process vs. Modern Solution

## Challenges

- Time consuming
- Unable to purify separate antigens
- Low purity
- Low yield
- Difficult to scale up
- Salt disposal issue

## Advantages

- Able to purify separate antigens
- High yield
- High purity
- Easy to scale up
- Time saving
- Environmentally friendly



# Process Highlights

1. Modern process to produce PT, FHA & PRN using bioprocess friendly, easily scalable, new generation chromatography platform.
2. Environmentally friendly.
3. Increase purity from 85% to >95%.
4. Reduce manufacture time from month to days.
5. Recovery increased from 10% to 30%.
6. Establish a sensitive, stable platform using Biacore to quantify PT & FHA.





# Modern Process for Meningococcal Vaccine



# Meningococcal Vaccine



13 clinically significant serotypes. A, B, C, W-135, Y responsible for 90% of global cases

Vaccine for A, C, W, Y are produced using capsular polysaccharide (PS), conjugant technology to enhance immunogenicity






# Traditional Process vs. Modern Solution

Traditional process challenge:


- Phenol use
- High-speed centrifuge



# Modern solution for Meningococcal Vaccine A,C,W,Y

Advantages vs. traditional process:

- No phenol use in process, benefit environment & operator's health & safety
- Easy to scale up
- Simple flow-through mode
- All 4 serotypes using same process
- Protein/DNA/endotoxin in products meet requirement



# Live Influenza virus production



# Influenza process overview



Scale-up from small scale to pilot scale in single-use format


Comparison of culture performance in 10 L and 50 L microcarrier culture in rocking bioreactors

Downstream purification in flow-through chromatography mode with Capto™ Q and Capto Core 700 chromatography media (resins)




# Cell growth in single-use bioreactor stage

## Cell concentration



## Average growth rate




Bead to bead transfer was successful and cell growth was comparable at 10 L and 50 L scale




# Virus growth kinetics

HA concentration and virus titer during culture




Cell morphology at time of harvest (96 h)



HA concentration at harvest was close to 12 μg/mL  
and the virus concentration was  $> 10^9$  infective units/mL




# Purification Workflow






# Purification results



Capto™ Q: Reduces host cell DNA




# Purification results



Capto™ Core 700: Reduces host cell protein



# Virus infectivity



Process does not impair virus infectivity



# Process summary

Estimation of doses per liter harvest, compared with WHO guidelines for protein and DNA impurities in influenza vaccine

	Split-inactivated vaccine <sup>1</sup>	Nasal LAIV <sup>2</sup>
Scale-up output/L harvest	175 doses á 15 µg HA	3075 doses á 10 <sup>7</sup> TCID <sub>50</sub> units
Harvest volume to produce 10 <sup>6</sup> doses	5760 L	325 L
Protein impurity <sup>3</sup>	30 µg protein/15 µg HA	1.5 µg protein/10 <sup>7</sup> TCID <sub>50</sub> units
DNA impurity <sup>4</sup>	3.0 ng/15 µg HA	0.15 ng/10 <sup>7</sup> TCID <sub>50</sub> units

<sup>1</sup> Split-inactivated vaccine contains 3 strains á 15 µg/HA (e.g., 3 × 15 = > 45 µg HA/dose á 0.5 mL).

<sup>2</sup> Comparison is based on a commercially available specification for a nasal LAIV. A dose of 0.2 mL contains 10<sup>7</sup> fluorescent focus units, which is assumed to be equal to TCID<sub>50</sub> titer.

<sup>3</sup> WHO guideline for protein impurity: max. 100 µg protein/strain

<sup>4</sup> WHO guideline for DNA impurity: < 10 ng DNA/dose = 3.3 ng DNA/15 µg HA.

Assuming a recovery of 25% for the overall process and a dose requirement of 10<sup>7</sup> TCID<sub>50</sub>, more than 1.5 million doses of monovalent live attenuated influenza vaccine could be produced from a 50 L cell culture



# Yellow fever virus propagation – from eggs to cells



# GMP manufacturing of viral vaccine



Xcellerex™ XDR-50 bioreactor

Vero cells (WHO-10-87)


- Cytodex™ 1 microcarrier
- Serum free, animal component-free medium

Yellow fever virus 17D



# Virus production drain down refeed

PFU equivalents from Eng and GMP bioreactor runs




The process consistency was high and virus titers were similar between runs





# Virus propagation and release of HCP

## HCP content after ELISA analysis




A feasible time for harvest is before the HCP peaks,  
to facilitate downstream processing



# Conclusions



# Conclusions



By modernizing legacy vaccine processes there can be improvements in:

- Yield
- Quality
- Scale-up
- Cost efficiency
- Environmental impact



# Acknowledgement

Dr. Tian, Yang Michelle (Fast Trak, GE Healthcare Life Science)

Dr. Zheng, Yongjian (Fast Trak, GE Healthcare Life Science)

Dr. Lu, Lifang Lillian (Fast Trak, GE Healthcare Life Science)

Tushar, Joglekar (Fast Trak India, GE Healthcare Life Science)

Chen, Yongshuo (GE Healthcare Life Science)

Liao, Yi Jazmine (GE Healthcare Life Science)

Li Xiaoling (Fast Trak, GE Healthcare Life Science)

Dr. Sui, Lili (Fast Trak, GE Healthcare Life Science)

Luo, Liang (GE Healthcare Life Science)

Xie, Hongyan (GE Healthcare Life Science)

Dr. Jia, Guodong (Fast Trak, GE Healthcare Life Science)

Yang, Hongyan (Fast Trak, GE Healthcare Life Science)

Dr. Lu, Lifang (Fast Trak, GE Healthcare Life Science)

Wang, Jiyuan (Fast Trak, GE Healthcare Life Science)

Yang, Jianjun (Fast Trak, GE Healthcare Life Science)

Du, Qizhen (Fast Trak, GE Healthcare Life Science)

Dr. Lundgren, Mats (GE Healthcare Life Science)

Dr. Berg, Mikael (GE Healthcare Life Science)

Larsson, Anders (GE Healthcare Life Science)



# Acknowledgement



## Uppsala team:

Lena Sandberg  
Peder Bergvall  
Camilla Estmer- Nilsson  
Christian Kaisermayer  
Ann-Christin Magnusson  
Christine Sund-Lundström  
Therese Lundberg  
Mats Lundgren  
Jakob Liderfelt  
Annika Morrison  
Elisabet Linde  
Åsa Eriksson  
Johanna Tschöp

## Marlborough team:

Geoff Hodge  
Patrick Guertin  
Ed Hayman  
Steve Turbayne  
Casey Cunningham  
Greg Mantenuto  
A. J. Theriault  
Chris Mach



# Thank you!

[www.gelifesciences.com](http://www.gelifesciences.com)

GE, imagination at work, GE monogram, and Global Star are trademarks of General Electric Company.

ÄKTA, Biacore, Capto, Cytodex, HyClone, ReadyCircuit, ReadyToProcess, Sephadex, ULTA , WAVE Bioreactor, and Xcellerex are trademarks of GE Healthcare Company or one of its subsidiaries.

HyperFlask is a trademark of Corning. Pluronic is a trademark of BASF.

© 2013-2015 General Electric Company – All rights reserved. First published Apr. 2013.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden



