

Introduction to process development for vaccine production

DCVMN 10 March 2017

Need for updated vaccine processing and process optimization for global access

Process development and optimization training workshop—DCVMN Taipei

Speakers



Dr. Mats Lundgren

Customer Applications Director, Life Sciences

GE Healthcare, Sweden



Mia Bennemo

Senior Research Engineer, Life Sciences

GE Healthcare, Sweden



29256323AA | Mar. 2017

Preliminary agenda

- Introduction to process development
 U for vaccine production
- Process economy
- Group discussion exercise—how to overcome technical and economical challenges in DCVMN companies
- Analytics
- Lunch break

- Upstream process development
- Downstream process development
- Quality by design (QbD) in process development
- Practical exercise: QbD
- Future scenarios, wrap up discussion and test



Introduction to process development for vaccine production



Outline

- Vaccine process history
- Why is process development for vaccines important?
- Vaccine processing
- Single-use technologies for vaccine manufacturing
- Conclusions









The evolution of vaccine processes

- First generation processes:
 - Focus on upstream, optional inactivation
- Second generation processes:
 - Separation based on centrifugation, filtration
- Currently developed processes:
 - Quality based approach: QbD
 - Focus on entire process including purification and virus safety



Historic vaccine timeline: propagation





Historic vaccine timeline: purification





The history of polio vaccines and GE

1955	1960	1960s	1970s	1980s	2010s
Inactivated polio vaccine (IPV) launched (Salk Type)	Attenuated polio vaccine launched (Sabin type)	Collaboration between Prof. Van Wezel (RIVM/NVI Netherlands) and GE (former Pharmacia) around microcarrier cultures of primary monkey cells	New IPV purification method using GE's chromatography resins	Switch to Vero cell production using Cytodex™ 1 microcarriers	Updating the IPV processes using modern technology from GE



Polio vaccine process



Production system	Vaccine type	Reference
Vero cell line	Polio vaccine (IPV). Naked virus (~ 30 nm). Type 1, 2, and 3 subtypes in the vaccine. (Sabine –OPV, Sabine-IPV, Salk-IPV)	Netherlands Vaccine Institute (NVI) - Vaccine 29, p.7188– 7196, 2011 - www.plosone.org, 1 December 2013, Vol. 8, Issue 12



Why is process development for vaccines important?



What are the challenges for vaccine producers?

Design of aged processes

- Many "weak steps", low yield, low robustness
- Lack of platforms, re-use of technology modules
- Open handling and regulatory concerns
- Regulatory practice does not support new technology implementation
- CAPEX demand very high due to weak processes
- Economy very dependent on scale

Adaptation to changing markets

- Markets for classic vaccines shrink in developed markets with high prices
- Need to remove hurdles for investment and improvement, including regulatory hurdles
- Reduce cost for highest standard production technology
- Overcome lack of flexibility in production infrastructure

Access to new vaccine technology

- Virus-like particles (VLP)
 - High safety
 - Low immunogenicity
 - Complex processes
- rec Antigens and adjuvants
 - Easy processing
 - Good safety
 - Immunogenicity dependent on adjuvant
- ...and more
 - Viral vectors
 - Plasmids, mRNA
 - Cells
 - Technologies in its infancy
 CAPEX = capital expenditure



29256323AA | Mar. 2017

What are the challenges for vaccine producers?





Process development—trends and solutions



- Quality must never be compromised
- You always need to spend enough time and money to understand what the process does to the product
- Your process will be most efficient in terms of time and cost, if you build its performance on solid understanding



Vaccines are difficult to characterize

Monoclonal antibodies	Vaccines	Implications
Often well-characterized	Often difficult to characterize	Less definitive analytical comparability pathways Less ability to monitor product quality in mid-process
Clear link to mechanism of action (MoA) and/or biomarker surrogate for clinical performance	Difficult to establish clinical potency surrogates	Challenging to improve process post-licensure
Consistent process and product	Sometimes more complex, less predictable process/product	Variability over product/process life cycle
Therapeutic patient population	Prophylactic patient population	"Process is product" philosophy to assure quality
Well-understood process; good detectability for test methods	Less understood process; difficult to measure attribute changes	Empirical process models for linking parameter inputs to quality outputs More stringent threshold for reporting manufacturing changes



A-VAX case study objectives I

Substantial changes in quality systems and regulatory approaches might be needed

- Apply QbD to develop a robust vaccine manufacturing process. This includes:
- Risk-based approaches to vaccine development
- Leveraging of science to gain process and product understanding
- Continuous improvement
- Merging of process and analytical controls for vaccine manufacturing
- Make the rationale for development more transparent in regulatory submissions



A-VAX case study objectives II

Substantial changes in quality systems and regulatory approaches might be needed

- Document techniques for safe and effective vaccines to reach market more quickly
- Strive to make reviews more efficient; decrease the number of postapproval supplements needed
- Develop realistic examples to better illustrate how QbD can be applied within the development space and overall product quality system
- Highlight and/or develop tools, frameworks, etc., to enable ICH Q8, Q9, Q10, and Q11 implementation strategies
- Tie key benefits with the strategies illustrated in the case study



Note of caution on A-VAX

It should be understood that that this document does not represent new regulatory policy, nor does it define a new "gold" standard for future regulatory submissions.

However, it is aligned with the available guidance from of ICH and other sources.

Individual companies will interpret and apply the principles differently. The extent of applicability will vary for each development effort. There are simplifying assumptions, e.g., the effect of multiple changes across unit operations is not considered.

There are aspects left out due to differences in opinion between participating companies. By no means should this case study be turned into regulatory expectations or standard.



Impact of process economy



Annual Report 2009: Crucell, incl. Berna Biotech

Typical vaccine business

 Legacy processes: manufacturing cost represents more than 50% of the revenue.

 Impact of process optimization on process economy (cost per dose reduction).



Vaccine processing



Vaccines and production





Diversification of technology—low efficiency





Vaccine production today

Processes developed
decades agoProcesses difficult to scale upOld cell substrates or
eggsCentrifugationLimited purificationFixed installationsSignificant expertiseRoller bottles

Unfavorable process economy		
Low yields		
Long process times		
Labor-intense processes		
Dedicated facilities		

Increased regulatory requirements Open handling

Batch variability

Serum supplementation



required

Vaccine production tomorrow

Processes developed decades ago	Processes difficult to scale up	Unfavorable process economy	
Platform cell lines	Scalable technologies enabled by, e.g., single-use technologies	Efficient and rational process design	
Efficient purification	<u> </u>		
		Elevible feeilitee	

Flexible facilites

Increased regulatory requirements

Closed handling

QbD

Chemically defined cell culture media



Vaccine manufacturing





Single-use systems suitable in vaccine production

- Vaccines often manufactured in relatively small batch sizes makes single-use technology appropriate
- Campaign manufacturing is common, single-use allows multiproduct manufacturing
- Pandemic preparedness requires faster development and manufacturing times
- Higher cost constraints on vaccine manufacturing call for improved process economics
- Safety concerns makes closed systems suitable



Single-use processing

Closed system processing—connecting upstream to downstream



- Standard or customized assemblies
- Considerations:
 - Sterility claims
 - Extractables/leachables
- Aseptic processing of large viruses (e.g., pox vectors)
- Improve economics
 - Reduced losses in sterile filtration

Quality must never be compromised

- You always need to spend enough time and money to understand what the process does to the product
- Your process will be most efficient in terms of time and cost if you build its performance on solid understanding



Need for updated vaccine processing and process optimization for global access.

gelifesciences.com

GE, the GE Monogram, Cytodex, FlexFactory, Sepharose, and Sephadex are trademarks of General Electric Company.

© 2017 General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

For local office contact information, visit **gelifesciences.com/contact**

GE Healthcare Bio-Sciences AB Björkgatan 30 751 84 Uppsala Sweden



