

**Millipore®**

Expert Pharm/BioPharm  
Products & CTDMO Services

# Empowering LMICs: Maximizing Impact through Collaborative Establishment of mRNA Platform Technologies

**DCVMN AGM**

Josephine Cheng  
Senior Modality Leader RNA Therapeutics  
Process Solutions, Merck Life Science  
Bali, Oct. 2025

**MERCK**



# Agenda

- 0 1 The world of collaboration
- 0 2 Highlighting 3 cases throughout the pandemic
- 0 3 Technical sharing of our recent publication
- 0 4 Summary



**Millipore®**

Expert Pharm/BioPharm  
Products & CTDMO Services

# Collaborations

01

**MERCK**



## Public-Private Partnerships (PPP)

COVAX Facility/AMC, pooled procurement, APAs

## Advance Purchase Agreements (APAs)

Upfront funding & dose commitments

## Pandemic Preparedness

CEPI & the 100 Days Mission; regional capacity

## Regulatory Pathways

WHO EUL (emergencies) & WHO  
Prequalification (PQ)

## Strategic Partnership

Acceleration of recourse allocation and maintain focus

## Technology Transfer (TT)

WHO/MPP mRNA CoE & partners; training & platform enablement

## CDMO

End-to-end dev & GMP manufacturing support

Millipore®



## Common Goals

vaccine collaborations aim to  
deliver faster, fairer, safer, and  
more sustainable access to  
immunization

MERCK



## Slide 4

---

**SL1**

What are these dots?

Sui-Ching Low, 2025-10-28T01:29:48.636



# How We Serve the Vaccine Industry

Millipore®


Inactivated/Live-attenuated Vaccines   Subunit Vaccines   Virus-Like-Particle (VLP) Vaccines   Viral Vector Vaccines   Plasmid DNA (pDNA) Vaccines   mRNA Vaccines   Bioteasing & CDMO   Emprove® Chemicals   Support

## Comprehensive Products & Services Portfolio

### Key Capabilities to Increase Efficiency and Safety of Vaccines

This document describes a range of products most commonly used for vaccine process development and manufacturing. Because each vaccine workflow is unique, however, selection of products most appropriate for your specific workflow requires testing and confirmation through MSAT support.

#### steps



#### example products



Millipore. SAFC. BioReliance.

Speed to clinic and simplified compliance

Collaborations  
Training



Industry and academic  
partnerships

MERCK



# Merck has long history of collaborating in vaccine industry

## Public Academia, industry references and collaborations

Millipore®

### Technology & services



### Platform development/ optimization



### Sponsorship & training





## Highlighting 3 cases throughout the pandemic

Rabies vaccine – Jenner Institute/ Oxford U. (2018)  
Shistosomiasis – Baylor College of Medicine (2019)  
WHO/MPP mRNA Tech Transfer Programme (2022)

02




# Collaborate to Build Effective and Cost-efficient Processes

## Early Phase Production for Ade-Vector-Based Vaccines

Millipore®

- **Objective:** Develop a cost effective, rapid platform for adenovirus vector using a Rabies vaccine candidate for clinical phase 1 material, using single use and filtration technologies.
- **Partner:** The Jenner Institute, University of Oxford, UK



Publication:

- Journal Paper: [Fedosyuk et al, Vaccine, 2019](#)
- Whitepaper: [jenner-vaccine-wp4773en-final.pdf](#)
- News: [Merck Supports Jenner Institute to Reach First Milestone in Covid-19 Vaccine Manufacturing](#)



### Key achievements

- Development of a platform at 4L batch scale for 2000 doses
- **1 week process (5 Days)**
- Optimized Clarification & TFF

14 Apr 2020

### Merck Supports Jenner Institute to Reach First Milestone in Covid-19 Vaccine Manufacturing

Merck and The Jenner Institute today announced that the Jenner Institute has laid the foundation for large-scale production of its Covid-19 vaccine candidate, ChAdOx1 nCoV-19.

were done hand in hand with Merck engineers



# Collaborate to Build Effective and Cost-efficient Processes

## Process optimization for rProtein-based vaccine candidates for neglected tropical disease Schistosomiasis



Millipore®


### Background:

- a disease caused by parasitic worms impacting 240M people.
- Since 2007, Merck works with WHO to provide up to 250 million praziquantel (PZQ) tablets annually, reaching 1.5B tablets to date.
- Global Schistosomiasis Alliance (GSA) coordinates efforts to eliminate schistosomiasis by advocating for resources and promoting effective prevention, treatment, and control strategies.
- The need for vaccines: targeting specific components of the schistosome parasites to elicit an immune response that can prevent infection, reduce need for praziquantel.
- Partner: Texas Children's Hospital Center for Vaccine Development (TCH-CVD), Baylor College of Medicine, USA

### Objectives:

Optimize a protein based vaccine Sm-TSP-2 derived from *P. Pastoris* for cost effective manufacturing targeting LMICs using Schistosomiasis candidate

Revised Clarification Process



### Key Achievements

- Reduction of process complexity
- Fouling & buffer use & membrane area requirements; Smaller manufacturing footprint
- 1.7 x higher concentration factor achieved; Shorter process time
- Replacing Imidazole by Histidine making it feasible for parenteral vaccine formulation.
- (2022) The work WAS leveraged for RBD-COVID-19 candidate development which was tech-transferred to Biological E (India) and Biofarma (Indonesia), and now approved to serve the world.**

Publication:

- Whitepaper Schistosomiasis process optimization: [neglected-tropical-diseases-wp12707en-mk.pdf](#)
- Whitepaper Baylor platform tech transferred for COVID vaccines: [vaccines-technology-transfer-wp12706en-mk.pdf](#)

**MERCK**

\*More details to check - <https://www.merckgroup.com/en/sustainability/health-for-all/schistosomiasis.html>




# Collaborate to empower LMICs for sustainable manufacturing WHO/MPP mRNA Tech Transfer Programme Programme

Millipore®

## Background:

The Programme is an initiative to facilitate the development and manufacturing of mRNA vaccines and therapeutics in low- and middle-income countries, launched in response to the COVID-19 pandemic to enhance global access by transferring mRNA technology and expertise.



MERCK




# ONE Merck global team supporting 15 partners Extensive Merck Support to the Programme and more to come



2022 2023 2024

2025





# Foundational Training Training

## Intensive hands-on training & openly discussion within the group where they all work on the same Programme

Millipore®



| Foundational Training Program                                |   |          |
|--|---|----------|
| Module   | Topic   | Duration |
| Module 1: Introduction to Millipore                          | Overview of Millipore and its products            | 15 min   |
|  | Introduction to Millipore's business model        | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's training program      | 15 min   |
| Module 2: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 3: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 4: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 5: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 6: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 7: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 8: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 9: Introduction to Millipore's products and services  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
| Module 10: Introduction to Millipore's products and services | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |
|  | Introduction to Millipore's products and services | 15 min   |

- Objectives: Skilled workforce enablement.
- 1<sup>st</sup> session Jan 6 – 17<sup>th</sup>, Singapore M Lab. 2<sup>nd</sup> session Oct. 27 – Nov. 7<sup>th</sup>, Molsheim, France M Lab.
- Participants from R&D, PD, Production and Quality department.
- High engagement and well appreciated by participants!
- Testimonial video by participants showing what they have learned in the training! <[LINK](#)>



# Visit us at the Booth! Recent publication related to mRNA



Millipore®

Whitepaper with Biovac: dedicated Biovac story on backwards integration strategy in the past two decade build up the Biovac today. (Jan, 2025)

Whitepaper with Afrigen: the journey to optimize GMP-ready TFF and Bioburden data for mRNA platform! (Aug. 2025)

Key mRNA technical references and key integrated solutions. (Oct. 2025)

Webinar in preparation (Dec. 11, 2025)




Merck




03

# Technical sharing on recent publication





The mRNA production process developed by Afrigen with various scales tested and support offered by Merck to optimize the process towards GMP-Ready process. **Millipore®**



Merck offers supports in (1) UF/DF and (2) bioburden control/ sterile filtration process optimization, (3) scale up optimization from 100ml to 1L, and (4) process conceptual design with single-use solutions.

(Abbreviations: IVT = in-vitro-transcription, BR=bioburden reduction, TFF= tangential flow filtration, UF= ultrafiltration, DF= diafiltration.)




## TFF operation options and parameters setting


TFF Filters with regenerated cellulose, 300KDa:

| IVT volume | Device mem. area   |
|------------|--------------------|
| 20 mL (PD) | 50 cm <sup>2</sup> |
| 100 mL     | 0.1 m <sup>2</sup> |
| 1 L        | 0.2 m <sup>2</sup> |

Permeate flux-controlled TFF:



Critical Flux Determination - 300KDa TFF Filter



- The permeate flux control can be achieved using a second pump or a regulation valve to maintain a consistent permeate flux, minimize fouling and concentration polarization, and is particularly helpful for sensitive or high-fouling applications.
- Initially both 100KDa and 300KDa membrane was tested, finally 300KDa was selected due to its superior flux.
- Flux excursion stopped at a flux of 80 LMH, operating flux range of 35-60 LMH was recommended.




# Scale up results for TFF1 to TFF4 from 100ml to 1L IVT scales demonstration robust scalability


- Results demonstrate good scalability of the TFF steps across TFF 1 to TFF 4, scaling from 100ml to 1L IVT Scale.
- Concentration – diafiltration – concentration, to optimize outcome
- Monitor TMP, fixed feed cross flow, and selected permeate flux from 35 – 60 for different TFF steps.

TFF1



100ml IVT Scale





1L IVT Scale





TFF2



TFF3



TFF4



Millipore®

MERCK



## Summary of TFF filtration results by average recovery and mRNA Integrity.

Millipore®

| Description           | 100 mL drug substance (n=3) | 1 L drug substance (n=3) |
|-----------------------|-----------------------------|--------------------------|
| <b>DRUG SUBSTANCE</b> |                             |                          |

### TFF filtration 1

|                           |                |                 |
|---------------------------|----------------|-----------------|
| % Recovery – TFF1         | 83.33 ± 2.31 * | 89.67 ± 4.51 ** |
| % Integrity – input TFF1  | 78.62 ± 2.10   | 78.43 ± 1.71    |
| % Integrity – output TFF1 | 76.63 ± 1.81   | 78.78 ± 2.07    |

### TFF filtration 2

|                           |              |               |
|---------------------------|--------------|---------------|
| % Recovery – TFF2         | 86.63 ± 2.36 | 95.68 ± 10.03 |
| % Integrity – input TFF2  | 73.27 ± 3.92 | 73.24 ± 2.99  |
| % Integrity – output TFF2 | 71.82 ± 0.66 | 70.03 ± 3.43  |

### TFF filtration 3

|                           |              |              |
|---------------------------|--------------|--------------|
| % Recovery – TFF3         | 91.53 ± 0.55 | 96.88 ± 3.01 |
| % Integrity – input TFF3  | 79.82 ± 1.35 | 73.31 ± 5.26 |
| % Integrity – output TFF3 | 76.07 ± 1.55 | 71.79 ± 2.91 |

| Description         | 100 mL drug product (n=3) | 3.7 L drug product*** (n=2) |
|---------------------|---------------------------|-----------------------------|
| <b>DRUG PRODUCT</b> |                           |                             |

### TFF filtration 4

|                           |              |                |
|---------------------------|--------------|----------------|
| % Recovery – TFF4         | 91.27 ± 4.97 | 100.25 ± 17.03 |
| % Integrity – input TFF4  | 76.27 ± 1.50 | 71.7 ± 1.27    |
| % Integrity – output TFF4 | 75.00 ± 2.07 | 68.40 ± 3.54   |



- Stable pre- and post-processing integrity results confirmed that the TFF process parameters did not shear the mRNA.
- In addition to the mRNA integrity not being affected by the TFF filtration, the process shows consistent high recoveries (>90% on average) for different scales and batches.
- These results demonstrate good scalability of the process and assure efficient concentration and buffer exchange.

MERCK



## Summary of Bioburden filtration results by average recovery and mRNA Integrity.

Millipore®

| Description                      | 100 mL DS (n=3)           | 1 L IVT scale (n=3) |
|----------------------------------|---------------------------|---------------------|
| <b>DRUG SUBSTANCE</b>            |                           |                     |
| Bioburden reduction filtration 1 |                           |                     |
| % Recovery – BR1                 | 98.21 ± 1.08              | 96.84 ± 8.18        |
| % Integrity – input BR1          | 76.63 ± 1.81              | 78.78 ± 2.07        |
| % Integrity – output BR1         | 76.63 ± 1.80              | 78.53 ± 3.06        |
| Bioburden reduction filtration 2 |                           |                     |
| % Recovery – BR2                 | 98.68 ± 0.79              | 86.45 ± 11.58       |
| % Integrity – input BR2          | 71.82 ± 0.66              | 70.03 ± 3.43        |
| % Integrity – output BR2         | 72.13 ± 1.03              | 71.02 ± 2.67        |
| Bioburden reduction filtration 3 |                           |                     |
| % Recovery – BR3                 | 98.82 ± 1.40              | 95.04 ± 5.05        |
| % Integrity – input BR3          | 76.07 ± 1.55              | 71.79 ± 2.91        |
| % Integrity – output BR3         | 76.42 ± 1.99              | 72.42 ± 1.12        |
| Description                      | 100 mL drug product (n=3) | 1 L scale * (n=2)   |
| <b>DRUG PRODUCT</b>              |                           |                     |
| Bioburden reduction filtration 4 |                           |                     |
| % Recovery – BR4                 | 98.57 ± 1.61              | 98.40 ± 0.37        |
| % Integrity – input BR4          | 75.00 ± 2.07              | 68.40 ± 3.54        |
| % Integrity – output BR4         | 76.40 ± 1.10              | 71.15 ± 2.90        |



- Recovery: Average recovery rates of above 95%, with minor losses attributed to filter wetting and flushing steps
- Integrity: mRNA integrity remained intact post filtration, with capillary gelelectrophoresis (CGE) confirming over 78% average integrity.
- Flux and pressure: At a constant flow rate of 300 LMH, the filter exhibited stable performance, handling feed volumes with minimal resistance.
- Based on the consistent performance of the bioburden reduction filtration results, the sterilizing-grade Millipore Express® SHC 0.5/0.2 µm filter can be suitably adopted for the sterile filtration step when needed.

MERCK



# Take Home messages



- 1 **Collaboration** formats can be **versatile**. Merck Life Science often collaborates through **scalable products** that include **quality dossiers, process optimization, scale-up implementation, technical training, and strategic partnerships**.
- 2 From the recombinant protein schistosomiasis vaccine and the viral vector platform rabies vaccine to COVID-19 mRNA vaccine candidates, **Merck has consistently been there** before and after the pandemic!
- 3 Merck supports the WHO/MPP mRNA Technology Transfer Programme with **global resources**, partners with the Center of Excellence **Afrigen** to optimize the platform, ensures that key purification steps are **GMP-ready processes**, and contributes to the goal of **empowering LMICs** to build a skilled workforce.



**Josephine Cheng**

Josephine.cheng@merckgroup.com

The vibrant M, Millipore, Mobius, Millistak, Durapore, Emprove, Pellicon, Eshmuno, Fratogel, Natrix, Millipak, Benzonase are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.


All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

**Acknowledgement**

Merck Team

Afrigen Team





## Vaccine ecosystem

### Collaboration type and key benefits

Millipore®

| Collaboration Type                 | Primary Goal   | Phase Focus                   | Key Benefits  | Typical Partners  | Examples                         |
|------------------------------------|--|-------------------------------|---|---|----------------------------------|
| Public-Private Partnership (PPP)   | Equitable access via pooled procurement & coordination | Access & delivery             | Portfolio aggregation, market shaping, equity         | Govts, WHO, Gavi, CEPI, UNICEF, industry, civil society | COVAX Facility & AMC             |
| Advance Purchase Agreements (APAs) | De-risk R&D/manufacturing; secure doses                | Pre-market financing & supply | Speed, risk-sharing, affordability                    | Governments/EC + manufacturers                          | EU APAs (e.g., AstraZeneca)      |
| Technology Transfer (WHO/MPP mRNA) | Build sustainable local mRNA manufacturing             | R&D → manufacturing           | Capacity building, autonomy, regional health security | WHO, MPP, Afrigen hub, LMIC manufacturers               | mRNA TT hub-and-spokes (Afrigen) |
| CDMO Partnerships                  | Provide development, scale-up & GMP manufacturing      | Development & manufacturing   | Speed to clinic/commercial, access to capabilities    | Sponsors + CDMOs  | End-to-end CDMO engagements      |
| Regulatory Pathways (WHO)          | Expedite emergency access; enable UN procurement       | Emergency & procurement       | Faster access, regulatory convergence                 | WHO, NRAs, sponsors                                     | EUL; Prequalification (PQ)       |
| Pandemic Preparedness (CEPI)       | Accelerate vaccine R&D for epidemic/pandemic threats   | R&D & response                | Funding, coordination, 100 Days Mission timelines     | CEPI, funders, R&D orgs, regulators                     | CEPI & 100 Days Mission          |
| Quality by Design (QbD)            | Design robust product/process quality                  | Development/CMC               | Reliable scale-up, regulatory-ready CMC               | Sponsor CMC teams; regulators                           | ICH Q8(R2)                       |
| Pharmacovigilance (PV)             | Ensure post-market safety & confidence                 | Post-market & programs        | AEFI detection, risk management                       | WHO, NRAs, UMC, health systems                          | WHO PV; VigiBase                 |

MERCK