
Tofflon Best Practice for Human Vaccine Manufacturing

www.tofflon.com









Tofflon Vision

What We Learn Today and Tomorrow in


The Science of Drug Manufacturing and Process Equipment Manufacturing

For the Future of


Smart Pharma Factory Building

Tofflon Innovation

Drug Development



Drug Manufacturing



Tofflon Innovative Solutions

Process Scale-up and Technology Transfer

Advanced Drug Manufacturing

Next Generation Drug Manufacturing

(Continuous Mfg. & Smart Mfg.)







Tofflon Innovations For Vaccine Manufacturing

www.tofflon.com

Tofflon Upstream Capabilities






Tofflon Downstream Capabilities





Tofflon Fill & Finish Capabilities







Washing & Drying

Rubber Cleaning

Filling

Capping


Inspection




Packaging www.tofflon.com

Integrated Solution









- Freeze Dryer
- Row-by-Row
- Isolator
- Auto Solution Preparation Systems (ASPS)
- Filling Line (Washer, Tunnel, Filler, Capper, External Washer)
- DCS



What Advantages KUFill Is Offering

- cGMP Compliance
- Modular Building Interfaces
- Automation and Isolation
- Systems Integration
- Modular Equipments
- Expandable Factory
- More New Technologies
- Integrated Project Management
- Added Value Technical Support and Service








Case study 1 : Recombinated Vaccine



Strengths :
Turnkey project
Integrated solution
Multiply planform
Disposable system
Fast Delivery
Flexible operation



Case study 2 : Inactivated Vaccine



Product : Inactivated Vius vaccine

Problem statement :

Closed and contained operation for handling of virus



- In production for virus inoculation, virus distribution, harvesting fluid handling, in QC for QC purpose viz virus sterility inspection, mycoplasma inspection, titer and other experiments
- Isolator designed with incubator cart to transfer material from isolator to incubator through RPT port.
- Considered sterility requirements and toxicity protection requirements for isolator.
- The internal dimensions of the isolator chamber are reasonably customized according to the different equipment, which are involved in different process.
- Separate air supply and separate air exhaust for isolator to ensure that isolator has no effect on the control of the differential pressure in the laboratory.

Tofflon Vision:

Smart Pharma Factory Builder

